

UniSart®
0.22 µm and 0.45 µm Nitrocellulose
Membranes for protein blotting



Sartorius Stedim Biotech at a glance

Sartorius was founded in Goettingen, in 1870 by Florenz Sartorius and started initially manufacturing short-beam analytical balances.

In 1929, Sartorius started the first ever industrial production of cellulose nitrate membranes. The Sartorius group has become an internationally leading laboratory and process technology provider with core areas of expertise in biotechnology and mechatronics.

Today, the recent combination of the Sartorius Biotechnology Division and Stedim, the pioneer in disposable aseptic bag systems, creates a global biopharmaceutical supplier. Sartorius Stedim Biotech is a leading provider of cutting-edge equipment and services for the development, quality assurance and production processes of the biopharmaceutical industry. Our integrated solutions covering fermentation, filtration, purification, fluid management and lab technologies are supporting the biopharmaceutical industry around the world to develop and produce

drugs safely, timely and economically. For next generation processes, we focus on single-use technologies and added-value services to meet upcoming requirements of the industry we serve. Strongly rooted in the scientific community and closely allied with customers and technology partners, we are dedicated to our philosophy of "Turning science into solutions".

Headquartered in Aubagne, France, and with main operations in Germany and the US, Sartorius Stedim Biotech owns production and sales companies worldwide and employs over 2,200 people. In the diagnostic field, based on its long experience of microporous membrane and its tradition of high quality, Sartorius Stedim Biotech has developed large pore size cellulose nitrate membrane, which are now a reference for consistency.





UniSart® 0.22 µm and 0.45 µm Nitrocellulose Membranes for Protein Blotting

UniSart® 0.22 µm and 0.45 µm Nitrocellulose Membranes for Protein Blotting

Line or dot blot strips are still widely used for in vitro diagnostic test.

Many infectious diseases generate a complex immuno response and necessitate checking the patient serum for the presence of several antibodies. For this purpose, the different Immunodominant antigen proteins have to be applied in line or dots onto a uniform binding substrate. Today, long ago after a first membrane was used in blotting, the nitrocellulose membrane is still the best substrate for proteins.

The UniSart® 0.22 µm and 0.45 µm nitrocellulose membranes have been designed to show exceptional consistency in IVD blot strips

All blotting applications

Gel blotting

The UniSart® 0.22 µm and 0.45 µm membranes are well adapted to all protein blotting systems such as; electro transfer, semi-dry or simply capillary blotting. It can be easily cut to the desired gel dimension. Under the electrical field and some capillary drive, proteins move from the gel to the UniSart® membrane which binds them immediately.

Line and dot printing

Today, most of the antigenic proteins can be obtained through recombinant techniques in purer state and larger quantities. As a result, these proteins can easily be dispensed in line or dots onto the membrane, rather than to be transferred from a gel. The smooth surface of the UniSart® 0.22 µm and 0.45 µm allows for an easy dispensing either with contact tip, in this case the hard and smooth UniSart® 0.22 µm and 0.45 µm surface is well adapted for that, either with non contact nozzles, allowing for sharp dot spots.

Consistency in your assay

Consistency

The request for consistency is the main one for IVD strip manufacturers. Only a very consistent nitrocellulose membrane, showing the same characteristics inside a batch and from batch to batch will allow the IVD manufacturer to maintain its production inside the targeted CVs.

The consistency of the final developed signal in a line or blot strip assay is closely dependent on the homogeneity of the membrane surface and Sartorius with its unique know how on membrane production has developed the UniSart 0.22 µm and 0.45 µm nitrocellulose membranes in order to show always the same characteristics inside a lot as well as inter lot. (see graphs)

A consistent membrane simply means:

- higher productivity
- reduced scrap rates
- reduced time to market
- reduced overall cost

Consistent wetting

Though nitrocellulose polymers are inherently hydrophobic, which contributes greatly to their very high protein binding capacity, the UniSart membrane has been specifically manufactured with hydrophilic reagents in order to show a fast and uniform wetting. The lines or dots that you will apply onto the membrane will always show the same shape. The best way to check it is to visually control the shape of the lines or dot immediately after dispensing. The wet lines or dots should then appear with clean and sharp edges. The width or size of lines and dots will be directly proportional to the dispensed volume.

Consistent binding

The UniSart 0.22 µm and 0.45 µm membranes shows an optimum structure for blotting. The small pore size allow for a huge protein binding capacity onto the sponge like 3 dimensional structure, up to 200 µg of equivalent IgG could be adsorbed in one cm² of UniSart membrane. At the same time the high ratio of pores allow for a fast flow through the membrane, as it can be noticed from the electron microscope photo.

The staining of the transferred proteins for control purposes can easily be done with standard methods like with Ponceau S, Coomassie Blue or Amido Black.


After transfer, the UniSart nitrocellulose can be blocked and rinsed according to the most usual protocol.

Consistent detection


After incubation with the test samples, the antibody-antigen complex can be revealed on the membrane strip following different methods. The most used being: colorimetric adding conjugate antibody with alkaline phosphatase or HRP for enzymatic detection, but also chemiluminescent and fluorescent with specific labeled molecules or even particles, as also used in lateral flow immuno assays.

Consistent sidedness

Due to the production process, the non supported nitrocellulose membranes like the UniSart 0.22 µm and 0.45 µm, always present two different sides. Both membrane sides should be first evaluated before selecting one.



Scan Electron Microscope photo of the Air side of the UniSart® 0.22 µm membrane × 4000 magnification



Scan Electron Microscope photo of the belt side of the UniSart® 0.22 µm membrane × 4000 magnification

Nitrocellulose production at Sartorius

Being the first ever manufacturer of cellulose nitrate membranes, Sartorius has many years of experience in the formulation of the lacquer as well as in the building of state-of-the-art casting equipment. In our new facility 2001 in Germany we have rigid environmental and process controls to provide a stable and clean manufacturing process.


Intra and inter lot consistency

The request for consistency is the main one for IVD strip manufacturers. Only a very consistent nitrocellulose membrane, showing the same characteristics inside a batch and from batch to batch will allow the IVD manufacturer to maintain its production inside the targeted CVs.

The consistency of the final developed signal in a line or blot strip assay is closely dependent on the homogeneity of the membrane surface and structure.

Sartorius with its unique know how on membrane production has developed the UniSart 0.22 μm and 0.45 μm nitrocellulose membranes in order to show always the same characteristics inside a lot as well as inter lot (see graphs).

Down and cross web thickness values of consecutive master rolls of UniSart 0.22 μm .



Graph Intra & inter lot consistency bubble point values on consecutive master rolls of UniSart® 0.22 μm membrane



of UniSart® 0.22 μm



Characteristics 0.22 μm and 0.45 μm :

The slightly more open structure will allow a faster flow through the membrane, with a limited decrease in total binding capacity. These two membranes 0.22 μm and 0.45 μm with slightly different characteristics will allow you to choose the best one for your particular assay.

Improved handling:

Both nitrocellulose membranes show an improved strength for ease of handling and cutting. They can also be easily glued with non solvent water resistant adhesives onto various plastic backing.

Versatile formats

Cut from wide master rolls, both membranes are available either in roll stock or as flat sheet. Any custom dimensions are available upon request.

Typical parameter value of the UniSart 0.22 μm and 0.45 μm :

	0.22 μm	0.45 μm
Flow rate for water:	27	70 [ml/(min.cm ² bar)]
Bubble point with water:	4.4	2.4 [bar]
Extractable content in water:	<1	<1 [%]
Thickness:	120	130 [μm]
Burst pressure:	0.8	0.2 [bar]
Wetting with water:	<1	<1 [secs]

Storage and handling of UniSart® membranes

Characteristics and certificate of quality

Perfect Quality assurance system is a prerequisite for manufacturing consistent membrane. Our QA encompasses all the membrane manufacturing steps. Detailed quality controls are performed from incoming materials, through casting process to the final release of a batch. Complete traceability of each UniSart membrane lot is then assured. The lot size is defined by the total membrane surface produced out of one mixing tank. All UniSart shipment is delivered with a certificate of quality.

Shelf life

The UniSart nitrocellulose membranes have a shelf life of 3 years after the Date of manufacturing. The end product manufacturer must determine the shelf life of the membrane once it is processed. Please refer to the Storage and Handling recommendations at www.sartorius-stedim.com/unisart.

Safety data

The unbacked nitrocellulose membrane material is considered for shipment as a dangerous good according to UN3270. Please refer to the Safety Data Sheet, available at www.sartorius-stedim.com/unisart.

Available materials and dimensions:

– UniSart 0.22 µm and 0.45 µm are available in roll stock or as flat sheet or disc.

Examples of packaging:

Product number	description
11327-160-110-N	UniSart 0.22 µm membrane sheets 160 mm × 110 mm box of 100
11327-----77	UniSart 0.22 µm membrane roll 300 mm on 50 m
11306-----77	UniSart 0.45 µm membrane roll 300 mm on 50 m



Sales and Service Contacts

For further contacts, visit www.sartorius.com

Europe

Germany

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen

Phone +49.551.308.0
Fax +49.551.308.3289

Sartorius Stedim Systems GmbH
Robert-Bosch-Strasse 5 – 7
34302 Guxhagen

Phone +49.5665.407.0
Fax +49.5665.407.2200

France

Sartorius Stedim FMT S.A.S.
ZI des Paluds
Avenue de Jouques – CS 91051
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.845619

Sartorius Stedim France SAS
ZI des Paluds
Avenue de Jouques – CS 71058
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.846545

Austria

Sartorius Stedim Austria GmbH
Franzosengraben 12
1030 Vienna

Phone +43.1.7965763.18
Fax +43.1.796576344

Belgium

Sartorius Stedim Belgium N.V.
Leuvensesteenweg, 248/B
1800 Vilvoorde

Phone +32.2.756.06.80
Fax +32.2.756.06.81

Hungary

Sartorius Stedim Hungária Kft.
Kagyló u. 5
2092 Budakeszi

Phone +36.23.457.227
Fax +36.23.457.147

Italy

Sartorius Stedim Italy S.p.A.
Via dell'Antella, 76/A
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41
Fax +39.055.63.40.526

Netherlands

Sartorius Netherlands B.V.

Phone +31.30.60.53.001
Fax +31.30.60.52.917

info.netherlands@sartorius.com

Poland

Sartorius Stedim Poland Sp. z o.o.
ul. Wrzesinska 70
62-025 Kostrzyn

Phone +48.61.647.38.40
Fax +48.61.879.25.04

Russian Federation

LLC "Sartorius ICR"
Uralskaya str. 4, Lit. B
199155, Saint-Petersburg

Phone +7.812.327.5.327
Fax +7.812.327.5.323

Spain

Sartorius Stedim Spain, S.A.U.
Avda. de la Industria, 32
Edificio PAYMA
28108 Alcobendas (Madrid)

Phone +34.902.110.935
Fax +34.91.358.96.23

Switzerland

Sartorius Stedim Switzerland AG
Ringstrasse 24 a
8317 Tagelswangen

Phone +41.52.354.36.36
Fax +41.52.354.36.46

U.K.

Sartorius Stedim UK Ltd.
Longmead Business Centre
Blenheim Road, Epsom
Surrey KT19 9 QQ

Phone +44.1372.737159
Fax +44.1372.726171

America

USA

Sartorius Stedim North America Inc.
5 Orville Drive, Suite 200
Bohemia, NY 11716

Toll-Free +1.800.368.7178
Fax +1.631.254.4253

Argentina

Sartorius Argentina S.A.
Int. A. Ávalos 4251
B1605ECS Munro
Buenos Aires

Phone +54.11.4721.0505
Fax +54.11.4762.2333

Brazil

Sartorius do Brasil Ltda
Avenida Senador Vergueiro 2962
São Bernardo do Campo
CEP 09600-000 - SP- Brasil

Phone +55.11.4362.8900
Fax + 55.11.4362.8901

Mexico

Sartorius de México S.A. de C.V.
Circuito Circunvalación Poniente
No. 149
Ciudad Satélite
53100, Estado de México
México

Phone +52.5555.62.1102
Fax +52.5555.62.2942

Asia | Pacific

Australia

Sartorius Stedim Australia Pty. Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175

Phone +61.3.8762.1800
Fax +61.3.8762.1828

China

Sartorius Stedim Biotech (Beijing) Co. Ltd.
No. 33 Yu'an Road
Airport Industrial Park Zone B
Shunyi District, Beijing 101300

Phone +86.10.80426516
Fax +86.10.80426580

Sartorius Stedim Biotech (Beijing) Co. Ltd.
Shanghai Branch Office
3rd Floor, North Wing, Tower 1
No. 4560 Jin Ke Road
Pudong District, Shanghai 201210

Phone +86.21.68782300
Fax +86.21.68782332 | 68782882

Sartorius Stedim Biotech (Beijing) Co. Ltd.
Guangzhou Representative Office
Unit K, Building 23
Huihua Commerce & Trade Building
No. 80 Xianlie Middle Road
Guangzhou 510070

Phone +86.20.37618687 | 37618651
Fax +86.20.37619051

India

Sartorius Stedim India Pvt. Ltd.
#69/2-69/3, NH 48, Jakkasandra
Nelamangala Tq
562 123 Bangalore, India

Phone +91.80.4350.5250
Fax +91.80.4350.5253

Japan

Sartorius Stedim Japan K.K.
4th Fl., Daiwa Shinagawa North Bldg.
8-11, Kita-Shinagawa 1-chome
Shinagawa-ku, Tokyo, 140-0001 Japan

Phone +81.3.4331.4300
Fax +81.3.4331.4301

Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.
Lot L3-E-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur, Malaysia

Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore

Sartorius Stedim Singapore Pte. Ltd.
1 Science Park Road,
The Capricorn, #05-08A,
Singapore Science Park II
Singapore 117528

Phone +65.6872.3966
Fax +65.6778.2494

South Korea

Sartorius Korea Biotech Co., Ltd.
8th Floor, Solid Space B/D,
PanGyoYeok-Ro 220, Bundang-Gu
SeongNam-Si, GyeongGi-Do, 463-400

Phone +82.31.622.5700
Fax +82.31.622.5799



◀ www.sartorius.com