



Sartopore® 2 0.1 µm

Sterilizing Grade and Mycoplasma Retentive MaxiCaps®, MidiCaps® and Capsules

MidiCaps® & MaxiCaps®



Description

Sartopore® 2 0.1 µm rated MaxiCaps®, MidiCaps® and Capsules are self-contained, ready-to-use membrane filter units for validated sterile filtration and reliable Mycoplasma removal in the pharma | biotech industry.

Applications

Typical applications include sterilizing grade filtration and Mycoplasma removal from:

- Animal Sera
- Cell Culture Media
- Media Components
- Bioprocessed Pharmaceuticals
- Prefiltration in front of virus filters
- Biological Fluids

and any other application requiring sub 0.2 µm filtration for enhanced sterility assurance.

Compatibility

The polyethersulfone membrane is compatible with a pH-range from pH 1 to pH 14 making Sartopore® 2 MaxiCaps®, MidiCaps® and Capsules ideal for filtration of solutions with high | low pH.

Easy to use

Sartopore® 2 MidiCaps® and Capsules are delivered as individually packed sterile units. On site, pre-use sterilization can be eliminated.

Flexibility

Sartopore® 2 0.1 µm MaxiCaps®, MidiCaps® and Capsules are available with various

filtration areas from 150 cm² | 0.16 ft² up to 1.8 m² | 19.4 ft² for easy adoption to any filtration process independent from the batch size.

Scalability

Consistent and predictable scale-up and down trials can reliably be performed as all Sartopore® 2 MaxiCaps®, MidiCaps® and Capsules are produced with the same type of membrane and materials and identical construction.

Cost Saving

The use of the disposable capsule design concept avoids investments into stainless steel filter housings and eliminates additional costs for cleaning of housings and cleaning validation.

Microbiological Retention

Sartopore® 2 0.1 µm MaxiCaps®, MidiCaps® and Capsules are validated as sterilizing grade filters according to ASTM F-838-05 standard and for Mycoplasma removal with a Log Reduction Value (LRV) of 7 for *Acholeplasma laidlawii*.

Quality Control

Each individual element is tested for integrity by Diffusion test prior to be released assuring absolute reliability.

Documentation

Sartopore® 2 MaxiCaps®, MidiCaps® and Capsules are designed, developed and manufactured in accordance with a ISO 9001 certified Quality Management System. A Validation Guide is available for compliance with regulatory requirements.

Specifications

Materials

Prefilter Membrane	Polyethersulfone, asymmetric
Endfilter Membrane	Polyethersulfone, asymmetric
Support Fleece	Polypropylene
Core	Polypropylene
End Caps	Polypropylene
Capsule Housing	Polypropylene
O-Rings	Silicone
Filling Bell	Polycarbonate

Pore Size Combination

0.2 µm + 0.1 µm

Available Sizes | Filtration Area

MaxiCaps®

Size 1	0.6 m ² 6.5 ft ²
Size 2	1.2 m ² 12.9 ft ²
Size 3	1.8 m ² 19.4 ft ²

MidiCaps®

Size 7	0.05 m ² 0.5 ft ²
Size 8	0.1 m ² 1.1 ft ²
Size 9	0.2 m ² 2.2 ft ²
Size 0	0.45 m ² 4.8 ft ²

Capsules

Size 4	0.015 m ² 0.16 ft ²
Size 5	0.03 m ² 0.32 ft ²

Available Connectors MaxiCaps®, MidiCaps®

S:	Tri-Clamp 50 mm (1 1/2")
O:	1/2" single stepped hose barb
F:	Tri-Clamp 25 mm (3/4")
H:	1/4" multiple stepped hose barb (only MidiCaps® Size 7 with filling bell)

Available Connectors Capsules

S:	Tri-Clamp 25 mm (3/4")
O:	1/4" multiple stepped hose barb

Operating Parameters

Max. Allowable Differential Pressure	4 bar 58 psi at 20°C (MaxiCaps®, Capsules)
	5 bar 72.5 psi at 20°C (MidiCaps®)
	3 bar 43.5 psi at 50°C (MaxiCaps®)
	2 bar 29 psi at 80°C (MidiCaps®, Capsules)
Max. Allowable Back Pressure	2 bar 29 psi at 20°C

Specifications

Extractables

Sartopore® 2 0.1 µm MaxiCaps®, MidiCaps® and Capsules meet, or exceed the requirements for WFI quality standards set by the current USP.

Regulatory Compliance

Individually integrity tested

Integrity test correlated to HIMA/ASTM F-838-05 Bacteria Challenge Test and Mycoplasma removal.

Non pyrogenic according to USP Bacterial Endotoxins

Pass USP Plastic Class VI Test

Non fiber releasing according to 21 CFR

Sterilization

Autoclaving:

134 °C, 2 bar, 30 min

No In-Line Steam Sterilization

Sterilization Cycles

Autoclaving: Min. 25

Technical References

Validation Guide

SPK5732-e


Order Information

Order Code	Pore size [µm]	Pack size [Pieces]	Test Pressure [bar psi]	Max. Diffusion [ml/min]
MaxiCaps®				
5441358K1--**	0.1	1	4.0 58	24
5441358K2--**	0.1	1	4.0 58	48
5441358K3--**	0.1	1	4.0 58	72
MidiCaps®				
5445358K7--**--A	0.1	4	4.0 58	4
5445358K8--**--A	0.1	4	4.0 58	6
5445358K9--**--A	0.1	4	4.0 58	9
5445358K0--**--V	0.1	2	4.0 58	18
Capsules				
5441358K4--**--B	0.1	5	4.0 58	1
5441358K5--**--B	0.1	5	4.0 58	2


** Connector Styles

(first code letter represents the inlet, the second the outlet)

Water Flow Rates for MidiCaps® with SS Connector



Water Flow Rates for MaxiCaps®



Standardized at 20 °C

Sartorius Stedim Biotech GmbH
 August-Spindler-Strasse 11
 37079 Goettingen, Germany
 Phone +49.551.308.0
 Fax +49.551.308.3289
 www.sartorius-stedim.com

USA Toll-Free +1.800.368.7178
 UK +44.1372.737159
 France +33.442.845600
 Italy +39.055.63.40.41
 Spain +34.90.2110935
 Russian Federation +7.812.327.5.327
 Japan +81.3.4331.4300

Specifications subject to change without notice. Printed and copyrighted by Sartorius Stedim Biotech GmbH. | W
 Publication No.: SPK2063-e13056
 Order No.: 85030-521-84
 Ver. 05 | 2013