
Installation Instructions | Installationsanleitung

Sartorius WZV 1.5 SE

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Special Information

These installation instructions describe the WZV1.5 SE series of weigh cells at the time of their latest revision, indicated by the revision number on the last page. However, the weigh cell(s) that you have received may differ. Any deviations are indicated in corresponding documents (quotations, performance specifications). Sartorius AG reserves the right to make changes to the technology, features, specifications and design of the equipment and the order data without notice.

The information included in these instructions was compiled with care. However, we cannot guarantee the applicability or correctness of this information. Because we are constantly developing and improving our equipment, it is possible that individual components or the instructions relating to them are not included in this manual in their entirety. It is also possible that some instructions in this manual do not relate to your model weigh cell.

If you believe that you are missing important information or have detected errors in the content of this manual, or if you would like more detailed information about individual components, please contact us:

Sartorius AG
Weender Landstrasse 94 – 108
37075 Goettingen, Germany
Phone: +49 (0)551 308-4440
Fax: +49 (0)551 308-4449
info@sartorius.com
www.sartorius.com

Intended Use

The WZV 1.5 SE weigh cell is a compact weighing system that can be installed in measuring and filling equipment that is used for small amounts of materials. Because of its compact design, the weigh cell allows several weighing procedures to be performed simultaneously within a small area.

Explanation of Symbols

- indicates handling instructions.
- indicates handling instructions that should only be carried out under certain conditions.
- > describes what happens once a handling instruction has been carried out.
- denotes an important point.
-  indicates a hazard.
-  indicates steps that must be performed with extreme care or those with special aspects that must be observed.

Definitions, Terms and Abbreviations

Calibration:

This describes the test to determine whether the output signal/readout on a weigh cell/balance is indicating the correct values. During calibration, a known reference weight is placed on the balance and the readout is compared to the actual weight of the reference.

Adjustment:

This describes the procedure to set the output signal/readout on a weigh cell/balance to the nominal output signal/nominal value of the readout/display. Using a known reference weight that has been placed on the balance, the (digital) setting of the measuring sensitivity is adjusted until the nominal weight is displayed.

EMC:

Electromagnetic force compensation

xBPI:

Extended Balance Processor Interface

Binary interface protocol used to control the functions of the weighing system.

PCB:

Printed circuit board

English – page 2

In cases involving questions of interpretation, the German-language version shall prevail.

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
Contents

2	Special Information	13	Communication
2	Intended Use	14	Regular Tests
2	Explanation of Symbols	14	Adjustment and Calibration
2	Definitions, Terms and Abbreviations	14	Calibration Procedure
4	Warnings and Safety Precautions	14	Reproducibility
4	General Installation Information	15	Adjustment Sequence
4	Sartorius Hotline	16	Calibration Sequence
5	Model Label	17	Specifications
5	Components	17	Mechanical Data
5	Weighing Unit	17	Metrological Data
5	Base	18	Electronics
5	Electronic Components	18	Electronic Connection
5	System Description	19	EMI / RFI Specification
5	Details	19	Changes and Modifications
5	Dimensions / Measurements	19	Service and Warranty
5	Weighing Pans	20	Dimensions (Scale Drawings)
5	Housing		
6	Base		
6	PCB		
7	Component Diagram		
7	Component Diagram for the WZV 1.5 SE Weigh Cell		
8	Installation		
8	Storage and Shipping Conditions		
8	Shipping, Damage Sustained During Shipment		
8	Incoming Inspection		
8	Equipment Supplied		
8	Handling Instructions		
8	General		
8	Cleaning		
9	Installing the Weigh Cell in a Housing		
10	Installation Diagram		
11	Operating the Weigh Cell		
11	Information on “Analytical Weighing” with Weigh Cells		
11	Handling Samples and Containers		
11	Weighing Electrostatically Charged Samples and Containers		
11	Weighing Magnetic or Magnetizable Samples		
12	Ambient Conditions		
12	Warmup Time		
12	Cooling		
12	Drafts		
12	Humidity		
12	Eliminating Electrostatic Charges		
12	Magnetism		
12	Oscillations		

Warnings and Safety Precautions

-  Always make sure the equipment is disconnected from power before performing any installation, maintenance or repair work.
 -  Only authorized service technicians who have been trained by Sartorius and who follow Sartorius' standard operating procedures for maintenance and repair work may open the equipment.
 -  Make absolutely sure to disconnect the electronics module from power before you connect or disconnect any electronic peripheral devices to or from the data interfaces.
 -  Do not use this equipment in hazardous areas where there is a risk of explosion.
 -  Do not use this equipment in areas where medical equipment is operated.
 -  Any incoming inspection or installation work that does not conform to the instructions in this manual will result in forfeiture of all claims under the manufacturer's warranty.
 -  Do not allow unauthorized or untrained personnel to handle or tamper with the weigh cell.
 -  If you use electrical equipment in installations and under ambient conditions subject to stricter safety standards, you must comply with the provisions as specified in the applicable regulations for installation in your country.
 -  Once it has been installed, the WZV 1.5 SE weigh cell must be checked for the following:
 - Compliance with directives and standards for electrical apparatus
 - Electromagnetic compatibility of the entire system*
 - Compliance with mandatory safety regulations
- * The term "system" as used in this manual refers to the entire array of installed components that is necessary to operate the WZV 1.5 SE weigh cell (housing, power supply, data interfaces, equipment used to transport items to and from the balance during weighing, etc.)

General Installation Information

- To prevent damage to the equipment, read these installation instructions thoroughly before using your weigh cell.
- These installation instructions describe only the general technological specifications of the weigh cell and the conditions that must be observed during installation.
- Do not expose the equipment to extreme temperatures, aggressive chemical vapors, moisture, shocks or vibrations.
- If there is visible damage to the equipment, disconnect the equipment from power.
-  The operator shall be responsible for any modifications to Sartorius equipment or connections of cables, or to equipment not supplied by Sartorius, and must check and, if necessary, correct these modifications (please refer to the section entitled "Service and Warranty" for more information).
-  On request, Sartorius will provide information on the minimum operating specifications (in accordance with the standards for defined immunity to interference).

If you have any technical questions regarding the construction, specifications or installation of the equipment, or if you encounter any defects in the weigh cell, please contact your supplier or the Sartorius hotline:

Weender Landstrasse 94-108
37075 Goettingen, Germany
Phone: +49 (0)551 308-4440
Fax: +49 (0)551 308-4449
info@sartorius.com
www.sartorius.com

System Description

Model Label

The weigh cells referred to in this manual are labeled with the model number WZV 1.5 SE.

Each weigh cell is clearly identified by a label indicating the model number (05.1).

The weigh cells are available in various versions, and the one that you have received may differ from this basic model number. If this is the case, a different model number will be indicated on the label and the weigh cell will usually feature a modified weighing pan designed for the associated applications. Any changes are indicated on the quotation or the performance specifications for the unit.

Components

Each weigh cell consists of two weighing modules, which can be operated independently of one another, and a built-in PCB.

The WZV 1.5 SE weigh cell consists of the following main groups of components (viewed from the top down):

Weighing Unit

The weighing unit is protected by a housing (03). It contains two mechanical weighing systems that work separately from one another. The corresponding weighing pans (01) and (02) are identified on the housing by numbers (03.4).

Base

The base (06) along with the insulating base (08) and the gasket (09) form the separation between the upper weighing unit and the electronic components underneath. In conjunction with the contact surfaces (06.1 and 06.2) and the compressors (07), the base (06) makes it possible to secure and position the weigh cell. The key edge (06.3) ensures the correct positioning of the weigh cell.

Electronic Components

The PCB (13) is held in place by a PCB frame (11) and protected by a partial cover (16). Power is supplied to the unit via the connector, which is also used for communication with the weigh cell.

Details

Dimensions / Measurements

Please refer to page 20 for the exterior dimensions of the weigh cell.

Weighing Pans

In the diagram of the weigh cell, the standard weighing pans (01) and (02) are shown. Sartorius adjusts the weighing pans at the factory using openings for this purpose (03.5.1-7 or 03.6.1-7). Once the pans have been adjusted, seals (04) are affixed over these openings.

If you require other weighing pans, Sartorius can help to develop and install these at the production facilities. Please observe the metrological requirements to ensure correct weighing results.

Any modifications made to the weighing pans without the prior consent of Sartorius or damage to the seals affixed over the adjustment openings will result in the forfeiture of all claims under the warranty.

The maximum overload in +Z direction (downwards) is 500 g. Overloads in this direction cause the weighing pans to sink by 0.5 mm (0.02") and touch the housing. The weighing units are not protected in -Z direction (upwards). The maximum allowable force is 1N (1 kg·m/s²).

Lateral forces (on the X- and Y-axis) acting upon the weighing pans are deflected.

Avoid inclines of over 1° in relation to the housing because this type of angle could damage the weigh cell beyond repair.

Housing

The housing (03) protects the weigh cell from damage caused by improper handling or dust. However, this housing cannot be compared to a conventional housing, which typically holds back much stronger forces from the sensitive weighing components inside. Despite the protective housing, you should handle the weigh cell carefully. Avoid subjecting the housing to strong forces, which could damage the weigh cell beyond repair.

The positioning pins (03.1 and 03.2) are used to adjust the weighing pan inside a complete housing that is not included in the equipment supplied. To make it easier to remove the weigh cell from the housing, notches have been milled into the positioning pins (03.1 and 03.2). A tool can be placed onto these notches to assist with pulling the weigh cell out of the housing.

The cone (03.3), which is placed at a slightly eccentric position, is also necessary for the positioning and orientation of the weigh cell in the housing.

The label with the model number (05.1) and the secondary number (05.2) are located on the left or right side of the housing.

Any damage to the warranty label (15) will result in forfeiture of all claims under the warranty.

Base

The base (06) secures the weigh cell to a mounting plate, which is not included in the equipment supplied. The contact surfaces (06.1 and 06.2) and the compressors (07) enable the weigh cell to be clearly positioned on the mounting plate. The key edge (06.3) prevents the weigh cell from being installed at a contorted angle on the mounting plate.

PCB

The PCB (13) with a partial cover (11 and 16) is connected to the base by an insulator (08). The insulator decreases the amount of heat transferred from the electronic components (dissipation loss) to the weighing components. A silicone gasket installed later on the insulating base prevents heat transfer through convection.

The PCB (13) contains all the circuits for both measuring channels and is affixed to the PCB frame (11) in a floating mount. The contact coil (14) provides the grounded connection between the two components.


The upper portion of the electronic components is protected by a cover (16), which is marked with an ESD label (17).

Commands are delivered and power is supplied to the unit via the plug connectors (13.2).

To make installation easier, the PCB frame (11) and PCB (13.1 and 13.2) are slanted. In addition, because the PCB is affixed in a floating mount, it is easier to install the weigh cell in a complete housing.

Component Diagram

Component Diagram for the WZV 1.5 SE Weigh Cell

Pos.	Component	Pos.	Component	Pos.	Component
01	Weighing pan 1	04	Seals	09	Gasket
02	Weighing pan 2	05.1	Model label	11	PCB frame
03	Housing / capsule	05.2	Secondary number	13	PCB
03.1	Positioning pin	06	Base	13.1	Tilted guide
03.2	Positioning pin	06.1	Contact surface	13.2	Plug connector
03.3	Cone for positioning	06.2	Contact surface	14	Contact coil
03.4	Weighing unit label	06.3	Key edge	15	Warranty label
03.5.1	Openings for adjustment	07	Compressor	16	Housing / cover
03.5.2		08	Insulating base	17	ESD label
03.5.3					
03.5.4					
03.5.5					

Installation

Storage and Shipping Conditions

 Once the equipment has been removed from the packaging, it may lose accuracy or become damaged if subjected to strong vibrations. Excessively strong vibrations may compromise the safety of the equipment.

- Do not expose the equipment unnecessarily to extreme temperatures, moisture, shocks or vibrations.
- Save the box and all parts of the packaging until you have successfully installed your equipment. Only the original packaging provides the best protection for shipment.
- Before packing your equipment, unplug all connected cables to prevent damage.
- Do not expose the equipment to acceleration in excess of 20 G ($20 \times 9.81 \text{ m/s}^2$).

Shipping; Damage Sustained During Shipment

The weigh cells are shipped in groups of 24 (one packaging unit contains 24 weigh cells). Repeat orders, including those for service purposes, can only be processed in the smallest unit for shipment (i.e., 24 weigh cells).

- The original packaging must be used to ship the weigh cells (including for return to Sartorius).

Incoming Inspection

- The OEM customer shall inspect the product and packaging immediately upon delivery for proper functioning, completeness, and absence of defects. This is to be performed in an incoming inspection within 10 days of delivery of the product. The incoming inspection must take place before the equipment is installed. Any obvious defects, errors or incorrect delivery must be reported in writing.

Equipment Supplied

There are no accessories included in the equipment supplied with the weigh cell.

Handling Instructions

General

 Since the weigh cells are not secured for shipment by transport locking devices, do not expose them to any conditions or forces in excess of those indicated in the section entitled “Mechanical Data” on page 17.

This is especially important when installing weigh cells in larger quantities or in filling equipment.

 In addition, do not subject the load receptors to forces in excess of those indicated under the heading “Mechanical Data” on page 17.

 The housing is made of thin aluminum and must be handled with care.

Before installing or removing a weigh cell, be sure to unplug the unit from AC power.

Be sure to take adequate grounding measures to protect the electronic components when handling the unit.

Cleaning

Please use care and observe the following rules when cleaning the weigh cell:

- Do not allow any substances or cleaning products to enter the weigh cell or come in contact with the electronic components of the unit.
- Observe the information included in the section entitled “Mechanical Data” on page 17.
- Use only cleaning products that will not damage the materials used in the construction of the weigh cell and the PCB/electronic components.
- Before cleaning, be sure to disconnect the unit from AC power.

Installing the Weigh Cell in a Housing

Before the functionality of the WZV 1.5 SE weigh cell can be put to use, a suitable housing must exist into which the weigh cell can be installed. This housing must have mechanical, electrical and electronic interfaces.

The design of the housing and the configuration of several weigh cells depend on the requirements of the OEM customer. Therefore, Sartorius does not offer this housing as part of the equipment supplied or as an accessory.

The following diagram (see page 10) of a configuration shows only a basic possibility of how WZV 1.5 SE weigh cells can be installed in a complete device. Detailed drawings necessary for the construction of this type of housing can be provided by Sartorius upon request. Sartorius can also act as a consultant for the construction of the housing so that disadvantages to and adverse effects upon the metrological qualities of the weigh cell can be avoided.

A basic design for the housing is shown on page 10. For simplification purposes, the enveloping, lateral housing is not pictured. On the diagram, two compartments can be clearly recognized: one for the electronic components (below; E5) and one for the weighing components (above; E2).

The weighing chamber, in which the weighing pans of the weigh cells are located (E1), is located above the weighing components.

Two plates separate both compartments. The upper plate (E3), preferably made of metal, holds and positions the weigh cell. The contact surfaces (06.1 and 06.2) and the compressors (07) hold the base (06) in an exactly defined position on the plate. There is an opening (E9) on the plate (E3) that is complementary to the corresponding key edge. This opening prevents the weigh cell from being installed in the wrong position.

The plate underneath, (E4), preferably made of plastic or another heat-insulating material, provides heat insulation between the lower and upper compartments.

At the same level as this plate are the insulating base (08) and a corresponding gasket (09), which shields the weighing components from a convection of warm air emanating from the lower part of the unit.

The backplane (below; E6) is the receptacle for the plug connectors (E7) as well as two springs (E8) for each weigh cell installed. The exact designation of the plug connector and the pin assignments can be found under the heading "Electronic Connection" on page 18. The two springs (E8) on the side are guides for the float-mounted electronic components.

The lower compartment (E5) (electronic components) can be actively cooled. However, be sure that cooling the compartment does not create high pressure or a vacuum, which would bypass the gasket (09) and lead to convection in the compartment above.


The compartment (E2) containing the weighing components has openings for the weighing pans and is open at the top. You should avoid actively cooling this area because convection along the openings for the weighing pans can lead to erroneous weighing results.

 Air flow in the area surrounding the weighing pan adversely affects the weighing results.

The upper plate facilitates the exact positioning of the weigh cell and several weigh cells in relation to one another and ensures a defined distance between the weighing pans. This positioning plate has a corresponding opening for the positioning pins (03.1 and 03.2) and a spring with an adapter, which presses on the cone (03.3).

 **Before you begin with the construction of this housing, please consult a Sartorius representative.**

Installation Diagram

Operating the Weigh Cell

Information on Analytical Weighing with Weigh Cells

Handling Samples and Containers

To avoid air buoyancy errors and fluctuations due to convection streams on the surface of the sample, the item(s) being weighed must be conditioned to the same temperature as the weigh cell.

These effects increase with the volume and surface of the sample. Therefore, it is important to choose a tare container that is the appropriate size for the sample.

Do not touch samples or containers used for samples with your hands. The hygroscopicity of any fingerprints that remain on the sample or container and temperature influences have an effect on the weighing results.

Carefully place the sample on the weighing pan, e.g., manually using forceps or automatically via a robot or a filling system.

When constructing a draft shield, you must take measures to keep temperature increases in the inner chamber of the unit to a minimum (e.g., bypass).

Weighing Electrostatically Charged Samples and Containers

When weighing electrostatically charged objects, substantial errors in measurement can occur.

This includes samples that display extremely poor conductivity (such as glass, plastic, or filters) because any charges that have built up – through friction, for example – can only be diverted through the weighing pan, and this takes a considerable amount of time.

The consequence is a force action between the charges on the sample and the fixed parts of the weigh cell or in the environment. This is noticeable in the constant change in the readout of the weigh cell.

In addition to purely mechanical measures (such as shielding the sample using a special weighing pan), neutralizing charges on the surface by bombarding them with ions of the opposite polarity is an especially efficient method for eliminating this effect.

The air surrounding the sample can become conductive through ionization. This allows charges to level out in air and dissipate with respect to the earth (ground).

Charges in the area surrounding the weigh cell can also cause considerable interference to the weighing procedure. If you construct a draft shield yourself, you must take appropriate measures against electrostatic charges.

Weighing Magnetic or Magnetizable Samples

For technical reasons, the use of magnetizable materials in the production of weigh cells cannot be avoided.

The functional principle of high-resolution weigh cells is, to a great extent, based on compensation of the load on the weighing pan by magnetic forces.

When weighing magnetic or magnetizable samples and containers, there can be an interaction of the parts of the weigh cell mentioned previously. This can cause the weighing result to be erroneous.

To reduce this effect, we recommend increasing the distance between the sample and the weighing pan using a non-magnetizable material (quadratic reduction of force with distance).

Magnetic interference to the weigh cells can also occur due to magnetic fields in the area surrounding the weigh cell.

Ambient Conditions

The table below shows the essential environmental parameters necessary to ensure correct weighing results.

Temperature range	20 – 25°C
Permissible temperature changes	+/-1 K/h
Moisture	45 – 50% rel. humidity
Vibration	Vibration should be avoided; vibration effects depend on the installation conditions and the filter settings
Draft	The system must be equipped with a draft shield
Electromagnetic fields	These must be avoided; please see pages 11 & 12
Electromagnetic charges	These must be avoided; please see pages 11 & 12
Magnetic fields	These must be avoided; please see page 12

Warmup Time

To deliver exact results, the weigh cell must warm up after it is connected to power for the first time.

The amount of warmup time required depends on the system in which the weigh cell is installed. This can only be determined after the weigh cell has been installed into the end device or system.

- ☞ Only after this time will the device have reached the operating temperature necessary for precise weighing.

Cooling

Differences in the ambient temperature can also lead to convection streams and should therefore be avoided.

- ☞ Heat generated by the electronic components must be dissipated via suitable cooling measures.
- ☞ The temperature difference between the mounting plate onto which the weigh cells are mounted and the ambient temperature may not exceed 0.5 K.

Drafts

The environment around the weigh cell must be protected from drafts.

- ☞ Drafts can lead to erroneous weighing results.

Humidity

The relative humidity may not exceed 45% – 50% (non-condensing).

- ☞ Fluctuations in humidity can lead to erroneous weighing results.

Eliminating Electrostatic Charges

Through appropriate measures, it is possible to ensure that the material being weighed, the weigh cell and all surrounding components have the same electrical potential (ground). The load receptor must rest on the grounded portion of the weigh cell.

- ☞ Electrostatic forces can lead to erroneous weighing results.

Magnetism

Make sure that there are no magnetic materials in the area around the weigh cell. Magnetic fields (especially alternating magnetic fields) can influence weighing procedures. Only non-magnetic adjustment or calibration weights may be used with the weigh cell.

- ☞ Magnetic forces can lead to erroneous weighing results.

Vibration

Do not expose the weigh cell to any vibration during weighing.

- ☞ Vibration can lead to erroneous weighing results.

Communication

The functions described in the following section entitled “Functions/Commands for the WZV1.5 SE Weigh Cell” are stored in the weigh cell electronics and can be used for communication with the weighing unit.

These functions are based on the Sartorius xBPI protocol. Upon request, Sartorius can provide a detailed general description of the xBPI protocol, which describes the entire range of functions. In the weigh cell, only the functions described in the following sections are implemented as a subset of the protocol. Sartorius warrants only the operation of the functions described here.

Communication is bidirectional via an RS485 interface without handshake signals.

On the plug connectors (13.2) there are interface connections for each weighing point on pins 1/2 or 11/12 (please see the section on page 18 entitled “Electronic Connections” for the pin assignments in the weigh cell PCB).

Connection to external devices (such as a PC) must be made “point-to-point.”

Functions / Commands for WZV 1.5 SE Weigh Cells

The recommended xBPI functions / commands are listed in the following table:

Functions	xBPI command	Input parameter	Output parameter
Read (net) value	0x1F	0	[Block of values]
Read S/N	0x100F	06	String5 (BCD-coded)
Weigh cell status	0x30	-	String8 (bit-coded)
Start calibration (not to start adjustment procedure)	0x28	0x75: Cal user weight 0x77: Cal factory weight	Status
Cancel calibration (cancels the setting)	0x29	-	Status
Reset a temp. error	0x41	-	Status
Reset	0x58	0x01	Status
Set UserCalWeight	0x100A	0x25 8 bytes, floating decimal (ANSI: double)	Status
Read UserCalWeight	0x1009	0x25	8 bytes, floating decimal (ANSI: double)
Toggle setting for “Ambient Conditions“ (for a definition of the filter settings, see (1), p. 69)	0x2C	0x01 0x02	
Zero the display	0x18		

Regular Tests

Repeated testing helps ensure correct weighing results are achieved.

We suggest the following test (see below).

This approach is based on a filling procedure in which a weighing container is used. Other applications should be approached accordingly.

- Read and store the value displayed while the weighing pan is empty. Calculate the difference between the values.
- Read the net value of the unfilled weighing container and store it. Then calculate the difference between the values.

The difference in the readouts must lie within a narrow range, which must be determined with intact cells installed in the complete system.

- Plausibility tests, for example, should be established for the measured weight and the duration of the filling procedure.

These regular tests set minimum requirements and must be performed to ensure that the weight values transferred from the weigh cell over the interface correspond to the actual weights being measured.

This type of testing also enables the user to better recognize systematic errors.

Sporadic Errors in Measurement

In general, there is no way to rule out the possibility of sporadic errors in measurement. Such errors also cannot be recognized through the tests described previously.

The possibility of this type of error occurring must be taken into consideration and evaluated in the overall context of the system and the process being performed. If applicable, this type of measuring error should be evaluated through a risk analysis of the entire system.

Adjustment and Calibration

Description:

Calibration and adjustment of the weigh cell must be performed separately for each weigh cell using an appropriate reference weight.

The weights used for calibration and adjustment:

- Must be within the calibration range (1g +/- 10 mg);
- Must have a high degree of accuracy so that they can be used for calibration;
- Must be able to be placed on the weigh cell in such a way that the unit achieves "stability," and the position of the calibration weight must correspond to the position of the actual material being weighed.
- The operator must observe any regulations regarding the traceability of the reference weights to a mass standard or other standards.
- It is especially important to implement measures to ensure that there are no magnetic or electrostatic interactions between the weights used and the weigh cell or filling equipment.

Calibration Procedure

- The calibration procedure in filling equipment must be carried out in accordance with Sartorius AG regulations.
- Calibration procedures are performed and monitored via the xBPI protocol.

Calibration should be done regularly (at least 1x a day) in conjunction with a reproducibility test.

If the weigh cell has not been operated for more than 1 hour (cold start), it needs to be calibrated. Before calibrating the unit, allow it to warm up (see page 12). If there are significant deviations, the weigh cell should be adjusted.

The frequency of calibration depends on:

- The operating requirements with regard to accuracy;
- The ambient conditions;
- The required degree of process reliability;
- Any other requirements as stipulated in the regulations.

Adjustment should be performed in the following circumstances:


- After the weigh cell has been installed in the system;
- Following repairs/maintenance on the filling equipment;
- If the filling equipment was out of service for a considerable amount of time;
- After calibration if it was determined that a weight fell outside of the tolerance range defined by the operator.

Reproducibility

The reproducibility test consists of a calibration procedure that must be performed repeatedly (10x) and the calculation of the standard deviation of the displayed net value, which is performed at the end.


 Only use adjusted weigh cells for weighing procedures.

Adjustment Sequence

	Overall status	Balance status	Balance status 2	Adjustment status
0x01	Weigh cell active	Invalid value		
0x02	Weigh cell in standby mode			Value too low
0x04	Weigh cell in Init mode			Value too high
0x08	Weigh cell adjusted	Stability		
0x10	Error			
0x20			Stability adjustment	
0x40		Command active		
0x80				

Calibration Sequence

Specifications

Mechanical Specifications

Dimensions		349.0 mm × 50.54 mm × 25.14 mm
Weight		approx. 500 g
Overload protection	in “-Z” direction in “+Z” direction	< 5 N < 1 N
Overload angle		< 1°
Positioning to the vertical axis		1/250 (e.g., 1 part in x and 250 parts) in Z-direction
Changes during measurement		< 1/10000
Max. acceleration of the WZV 1.5 SE weigh cell		20 G (=20 × 9.81 m/s ²)

Metrological Specifications

Capacity	1 g
Max. readout	1.02 g
Interface resolution	0.01 mg
Reproducibility	≤ 0.05 mg
Linearity (Lin)	0.1 mg
Temperature coefficient of sensitivity	20 ppm/K (within the indicated temperature range)
Weight of the weighing pan	5 g ± 10 mg
Required calibration weight	Nominal weight: 1 g ± 10 mg

(Measured as standard deviation: 10x with a 1g weight placed in the center of the weighing pan; no outside influences, e.g., draft, electrostatic charges, etc.)

Two weights of 500 mg each are placed on the weigh cell.

- a: Place weight #1 on the weigh cell: read value 1
- b: Remove weight from weigh cell
- c: Place weight #2 on the weigh cell: tare the unit
- d: Add weight #1 to the weigh cell in addition: read value 2
- e: Calculate the linearity using $\text{LinError} = 1/2 \cdot \text{abs}(\text{value 2} - \text{value 1})$

The calibration uncertainty of the calibration weight may be 1/3 of the allowable weighing uncertainty at 1g. The properties must be comparable to class E2 weights according to OIML.

Electronics

Data interface	Serial RS-285 interface
Connection	Point-to-point connection
Protocol	Sartorius xBPI protocol ¹⁾ ; asynchronous semi-duplex
Handshake	No
Character coding	8-bit binary data
Parity	Odd
Synchronization	1 start bit, 1 stop bit
Check sum	See[1] ²⁾ , page 12
Transmission rate	25 measured values per second
Transmission speed	38,400 baud
Application programs	None
Requirement for external power supply	+9 V +15..30 mA +9 V +15..30 mA +6 V +60..90 mA
Ripple	max. 200 mV for these voltages

¹⁾ "Description of the Protocols and Functions of the Weighing Platform Interface (xBPI Description)"

²⁾ When switching on the supply voltage, be sure to follow the correct order. First, switch on the voltage for -9 V and +9 V analog components in a synchronized manner. If necessary, switch on +9 V first and then -9 V. Afterwards, switch on the voltage for +6 Volt digital components. During this phase of switching on the voltage, currents of up to 100 mA can occur.

Electronic Connection

The supply voltage and the interface are connected via the contact surfaces of the PCB (13.2) (see page 7).

The backplane, which is not included in the equipment supplied, holds the PCB plug connectors (Syscon TEK Corp., model 1101)^(*) used to connect the contacts).

Pin Assignment Chart of the Weigh Cell PCB: (Pin 1 is marked)

Solder side of the PCB	Pin			Pin	Components of the PCB
TxD_P_CH1	2	0	0	1	TxD_N_CH1
"Free"/not for ext. use	4	0	0	3	"Free"/not for ext. use
+9 V	6	0	0	5	GND
+6 V	8	0	0	7	-9 V
"Free"/not for ext. use	10	0	0	9	Master Reset MR
TxD_N_CH1	12	0	0	11	TxD_N_CH1

^{*)} Edge card connector 1101 S series

Sartorius WZV 1.5 SE weigh cells are components that are intended for installation in a complete system.

Therefore, only very basic measures have been implemented against electromagnetic/high-frequency interference.

Comprehensive protection against electromagnetic interference should be provided for the complete system into which the weigh cells are to be installed.

Appropriate measures must be taken to prevent interference that is especially conduction-bound.

A complete test of the functions of the weigh cell can only be performed once the weigh cell has been installed in a system.

Changes and Modifications

Sartorius AG reserves the right to make changes to the technology, functions, and design of the instrumentation without notice, unless otherwise specified in this document.

Unless otherwise agreed upon in writing, equipment from Sartorius AG features a warranty in accordance with the valid "General Terms and Conditions" and the Terms and Conditions for OEM products provided by Sartorius AG.

This warranty covers malfunctions and defects caused by failures in the construction, workmanship or materials used in production.

The warranty does not cover:

- Interference or damage resulting from improper use or mishandling of the equipment,
- Parts that have become defective as a result of normal wear and tear,
- Consumables.

All claims under the warranty will be forfeited if:

- The user or any unauthorized person makes technical changes to the equipment,
- Any equipment that has not been authorized by Sartorius AG for use with WZV 1.5 SE weigh cells is used in conjunction with the unit, and damages result in the weigh cell or the system into which it has been installed.

Returns

Returned weigh cells must be clean, in a completely hygienic state and carefully packed in the original packaging. Any parts that were contaminated by media or media components must be cleaned, disinfected and sterilized in accordance with the safety regulations that are valid for the particular application area, e.g., for chemical or biological safety.

The sender must provide evidence that safety regulations were adhered to while the weigh cell was used. For this purpose, a signed decontamination declaration must be included with the weigh cell being returned. Without proof of decontamination or a description of measures that were performed to ensure safety, weigh cells will not be repaired.

The sender will be responsible for any damage sustained during transport and any necessary cleaning or disinfection that must be done once the weigh cell has been received.

Within the scope of product service, only complete weigh cells can be exchanged.


Maintenance and repair work can be performed only on complete weigh cells at Sartorius service centers.

If you require service, please contact your authorized service representative or Sartorius AG directly at the address below:

Sartorius AG
Northern Service Center
Weender Landstrasse 94-108
D-37075 Goettingen, Germany
Phone: +49 (0)551 308-3729 / 3740 / 3741
Fax: +49 (0)551 308-3730

Dimensions (Scale Drawings)

Dimensions of the WZV 1.5 SE weigh cell

☞ Supplemental documentation is available from Sartorius AG upon request.

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Sartorius AG
Weender Landstrasse 94-108
37075 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius.com

Copyright by Sartorius AG,
Goettingen, Germany.

All rights reserved. No part of
this publication may be reprinted
or translated in any form or by any
means without the prior written
permission of Sartorius AG.

The status of the information,
specifications and illustrations in
this manual is indicated by the
date given below.

Sartorius AG reserves the right to make
changes to the technology, features,
specifications and design of the
equipment without notice.

Status:
November 2005, Sartorius AG,
Goettingen, Germany

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☞ Ergänzende Unterlagen können bei der Sartorius AG angefordert werden.