User Manual | Benutzerhandbuch | Mode d'emploi

Benchtop Freeze-Thaw System

Table of Contents

English page 2

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Introduction	3	CryoPilot Software	11
CryoVessel	3	Overview	11
Celsius®	3	System Requirements	11
References	3	Components	12
Safety	3	Installation	12
General Precautions	3	User Interface	13
Heat Transfer Fluid	4	Getting Started	17
CryoPilot A	4	Program Flow	17
CryoMixer Jr.	4	Reviewing Data	17
		Using the Chart	18
System Components	4	Exporting Data	18
CryoPilot A	4	Printing	18
Computer	4	Profiles	18
CryoWedge Module	5	Data Files	20
CryoCassette	5	Diagnostics Troubleshooting	21
S ³ Module	6	21 CFR Part 11	22
Celsius®-Paks	6	System Operation	23
CryoMixer Jr.	6	Preparing for a Run	23
CryoHose	6	Running a Profile	25
CryoWrap	6	Removing Frozen Celsius®-Paks from the S ³ Module	25
Data Acquisition Device	6	Dispensing Thawed Samples from the S ³ Module	25
System Installation	7	System Maintenance	26
Tools Equipment Required for Installation	7	Cleaning	26
System Component Inventory	7	Changing the Heat Transfer Fluid	26
Selecting a Location CryoWedge and CryoHose Installation	7	Appendix	27
S ³ Module and CryoHose Installation	8	Information and Instructions on Disposal and Repairs	27
CryoPilot A Setup	9	Contact Information	27
Heat Transfer Fluid Fill	9		
CryoMixer Jr. Setup	10		
Computer Setup	10		
Data Acquisition Setup	11		

Introduction

Sartorius Stedim Biotech offers scalable freeze-thaw solutions using patented technologies for the handling, storage, and transport of biopharmaceuticals, including tools for process development, pilot scale operation, and commercial scale production. The benchtop system configurations are designed to reproduce the controlled freeze-thaw processes that occur in the large-volume, production-scale Celsius® and CryoVessel processes. CryoVessel is a stainless steel product line, and Celsius® is a single-use product line. Both systems use the CryoPilot A chiller and CryoMixer Jr. shaker.

CryoVessel

The bench-scale system for CryoVessel is called the CryoWedge. The CryoWedge module duplicates a small section of the CryoVessel by matching the heat exchange surface geometry, freeze path length, and material of construction.

Product volume requirements for feasibility studies can be further reduced by the use of the CryoCassette. These are 30 mL and 100 mL stainless steel containers that can be fitted into the CryoWedge in polymeric sleeves.

Celsius®

The bench-scale system for Celsius® is called the S³ Module; S³ stands for Scale-down, and Stability Studies. The S³ Module models the FT16 and FT100 Freeze-Thaw Modules by matching the freeze path length. The Celsius® S³ Freeze-Thaw system produces similar freeze-thaw times allowing for scalable freeze-thaw tests and stability studies.

References

Refer to the following documents for more information about the components of your benchtop freeze-thaw system:

- Huber Unistat 705 Operating Instructions
- Huber Unistat 705 Electrical Schematics
- Dow Syltherm HF Product Information
- Dow Syltherm HF Material Safety Data Sheet (MSDS)
- DuPont Suva 507 Material Safety Data Sheet (MSDS)
- DuPont Suva 23 Material Safety Data Sheet (MSDS)

Safety

Carefully follow all instructions to prevent serious injury or death to personnel and damage to equipment. Any use of the system components in applications for which they were not intended voids the units' warranties and releases Sartorius Stedim Biotech from any and all liability.

The following symbols indicate important information:

Warning:

Potentially hazardous situations that, if not avoided, could result in death or serious injury.

Caution:

Potential for minor personnel safety hazards or of potential damage to equipment or property.

Note:

Information that is important but not safety-related.

General Precautions

Ensure that company safety policies are followed at all times.

Disconnect power and follow proper lock-out | tag-out procedures when performing maintenance inside system components.

Do not stretch the CryoHoses in order to attach them to the chiller or the application. Instead, move the chiller or application in order to arrange the hoses in a relaxed position. Allow adequate slack to prevent kinks.

System components are heavy. Follow proper lifting and rigging techniques.

System Components

Heat Transfer Fluid

The system uses Dow Syltherm HF silicone heat transfer fluid (HTF). Do not use any other heat transfer fluid.

There are several hazards associated with use of HTF:

- Temperatures as low as -75°C and as high as 250°C can cause severe injury to skin, eyes, and other soft tissue. Use appropriate personal protective equipment when working with HTF or near exposed hot or cold equipment surfaces.
- Syltherm HF has a closed-cup flash point of 63°C (145°F) per ASTM D92. According to NFPA and OSHA, it is classified as a Class IIIA combustible liquid. Follow proper handling and storage procedures.
- Syltherm HF has a low surface tension and it can be difficult to avoid spills when filling or draining the system. Once spilled it is extremely slippery. Use absorbent material such as polypropylene pads to keep the bench and floor free of spilled HTF.
- Waste HTF and absorbent must be disposed of as hazardous waste. Refer to the product information and MSDS for more information.

CryoPilot A

Only qualified personnel may service the refrigeration equipment. The CryoPilot A uses R-507 and R-23, which are hydrofluorocarbon (HFC) refrigerants. These refrigerants are subject to environmental regulations and must not be vented to the atmosphere. Refer to the MSDS for more information.

CryoMixer Jr.

During operation, keep hands clear of the top cover plate. Do not exceed the maximum recommended agitation speed listed in the System Installation section.

CryoPilot A

The CryoPilot A is a recirculating chiller that controls the temperature of a heat transfer fluid (HTF) to a temperature setpoint.

A cascade refrigeration system cools the HTF when its actual temperature is above the temperature setpoint. The main components of the refrigeration system are a high temperature compressor, a low temperature compressor, a condenser, an automated expansion valve and an evaporator. The refrigerants used in the high and low stages are R23 and R507, respectively. The refrigerant condenser is air-cooled.

A pair of electrical resistance heaters heats the HTF when its actual temperature is below the temperature setpoint. A 1500W large heater is used for large temperature changes and a 100W fine heater is used for small temperature changes. This dual heater arrangement provides for precise control of HTF supply temperature.

A variable speed pump circulates the HTF through the evaporator, the heaters and the external application. For CryoWedge applications, the HTF flows through the jacket and the core heat exchanger area. For S³ Module applications, the HTF flows through the walls adjacent to the Celsius®-Pak clamps. An expansion reservoir simplifies priming the system and allows for HTF expansion and contraction.

The CryoPilot A is microprocessor controlled. The controller user interface has several pushbuttons, a rotary encoder, an LCD display with touch screen, and a system diagram with status LEDs.

The following are key specifications for the CryoPilot A:

- Electrical Power Supply
Note: The CryoPilot A is delivered without a plug. Due to the variation in local electrical receptacles, the customer must supply the plug.
 - North America:
208 VAC, 60 Hz, 1~, 20 A
 - Europe:
230 VAC, 50 Hz, 1~, 15.6 A
- Fluid Temperature Range:
 -75°C to 250°C
- Dimensions W x D x H:
 $425 \times 400 \times 720$ mm
($16.7'' \times 15.7'' \times 28.3''$)
- Weight: 90 kg (198 lbs)

Computer

The CryoPilot A is controlled and monitored using CryoPilot software on a computer. The CryoPilot software controls and monitors the benchtop system components including chiller, mixer, and thermocouples. See the section entitled CryoPilot Software for minimum system requirements and additional information.

Computers purchased from Sartorius Stedim Biotech as part of a benchtop system have the following specifications:

- Laptop form factor
- Operating System:
Microsoft Windows 7 Professional, 32-bit
- Microsoft Word and Excel
- Adobe Acrobat
- CryoPilot software pre-installed

CryoWedge Module

CryoWedge module names correspond to the volume of production scale CryoVessels. For example, a CryoWedge 125 (20") is meant to simulate a 125L CryoVessel.

- Working Volume: 350 – 4000 mL, depending on the size of the CryoWedge.
- Materials of Construction:
Body: 316L stainless steel
Gaskets: silicone
Cover: transparent polycarbonate
- Connections: HTF in, HTF out with quick-disconnect couplings
- Weight, Empty:
 - CryoWedge 20 (12"):
15 kg (32 lbs)
 - CryoWedge 125 (20"):
21 kg (46 lbs)
 - CryoWedge 300 (30"):
44 kg (97 lbs)
 - CryoWedge 300 (34"):
49 kg (108 lbs)
- Size Equivalence:
 - CryoWedge 20 (12"):
CryoVessel 20 L – 40 L
 - CryoWedge 125 (20"):
CryoVessel 60 L – 150 L
 - CryoWedge 300 (30"):
CryoVessel 200 L – 300 L
 - CryoWedge 300 (34"):
CryoVessel 300 L

CryoCassette

The CryoCassette minimizes the amount of product used within the CryoWedge. The CryoWedge can then be filled with buffer while the CryoCassette is filled with product. The design of the CryoCassette maintains the freezing and thawing characteristics of the CryoWedge to ensure scalability to the production CryoVessels. The CryoCassette also serves as a scalable stability container.

Figure 1:
Stainless Steel CryoCassettes

- Working Volume:
10 – 100 mL, depending on the size of the cassette.
- Available in 30mL and 100mL sizes.
- Materials of Construction:
 - Body: 316L stainless steel
 - Boat: polypropylene (30 mL) or polyethylene (100 mL).
- Sterilization:
Customer autoclavable.
- Finish: 10µ-in Ra, EP

S³ Module

The S³ Module freeze-Thaw Module provides the capability to model Celsius® FT100 and FT16 systems with as little as 20 mL of product. The configuration of the heat transfer surfaces within the module provide the same bi-directional heat flux as seen in the FT100 system utilizing frame and Celsius®-Pak designs.

Figure 2: S³ Freeze-Thaw Module

- Capacity: 1 – 10 product samples
- Working volume: 20 mL to 1 L, depending on size of Celsius®-Paks
- Materials of construction:
 - Shell: 302-304 stainless steel Shell
 - Base: 6061-T6 aluminum
 - Shell frame cover:
5052 aluminum
 - End plates: Rigid PVC foam sheet
 - Internal shell insulation:
Silicone foam rubber
 - HTF connection end cover:
Delrin
 - Cover: Transparent PVC
- HTF piping maximum design pressure: 4.1 barg (60 psig)
- Connections:
HTF in, HTF out with installed quick-disconnect couplings; drain
- Weight:
Approximately 18 kg (40 lbs)

Celsius®-Paks

All Celsius®-Paks, regardless of size, have the same material of construction and the same freeze and thaw path length as shown in Figure 3. Since Celsius®-Paks are single-use, the need for cleaning is eliminated.

Figure 3:
Comparison Between 30 mL and
16.6 L Celsius®-Paks

- Available in 30 mL and 100 mL sizes
- Ready to use, single-use collapsible container
- Pre-sterilized (25–45 kGy gamma radiation)
- RF sealing compatible
- Designed to allow for sterile temperature measurement
- Celsius®-Paks overwrap pouches are available. The use of an overwrap is recommended when conducting studies at higher temperatures (room temperature and above).
- 30 mL and 100 mL Celsius®-Pak Materials of construction:
 - Product contact layer: EVAM (ethylene vinyl acetate copolymer, monomaterial)
 - Gas and moisture barrier layer: EVA/EVOH/EVA
 - External robust, handling layer: EVA
 - Fill and drain ports: EVA
 - Extension tubing: C-Flex®
 - Luerlock®: Polycarbonate or Polypropylene (C-Flex® option)
 - Double barb elbow or union: Polypropylene (C-Flex® option)
 - Thermowell: Pt-cured silicone

CryoMixer Jr.

The CryoMixer Jr. is a reciprocating mixer that is positioned below the CryoWedge or S³ Module in order to ensure homogeneity of the thawed product. The reciprocating motion starts automatically when activated by the CryoPilot software. The CryoMixer Jr. is used for both CryoWedge and S³ Module agitation applications.

Figure 4: CryoMixer Jr.

Specifications:

- Stroke: 1" (25mm), reciprocating
- Electrical Power Supply:
 - North America: 115 V, 60 Hz, 1~
 - Europe: 230 V, 50 Hz, 1~
- Dimensions W x D x H: 505 x 585 x 120 mm (20" x 23" x 5")
- Weight: 24 kg (52 lbs)

CryoHose

All systems use two CryoHoses for the supply and return of HTF from the CryoPilot A unit. These hoses are stainless steel with surrounding insulation.

CryoWrap

CryoWrap is a neoprene sleeve that, when wrapped around CryoHose connections, prevents frost build up during CryoWedge or S³ module freezing.

Data Acquisition Device

The data acquisition device is compact module used to digitize thermocouple signals for monitoring by the computer. The data acquisition device has a breakout cable with receptacles for up to 16 T-type thermocouples.

System Installation

This chapter will guide you through the setup of the benchtop system.

Tools | Equipment Required for Installation

The following tools are needed to set up the Benchtop Freeze-Thaw System:

- Assortment of cross head and flat head screw drivers
- Vise grip wrench with a minimum 6" (15 cm) handle
- Assortment of box or crescent wrenches
- Absorbent material such as polypropylene pads for cleaning up any spilled HTF
- Multimeter suitable for measurement of supply voltage

System Component Inventory

When all the equipment containers are unpacked, confirm that the following parts are available for installation:

- CryoPilot A
- CryoWedge or S³ Module
- CryoMixer Jr.
- (2) CryoHoses
- (4) CryoWraps
- T-Type thermocouples (up to 10 depending on system)
- USB to CryoMixer Jr. cable
- Computer with CryoPilot Software
- National Instruments USB-9213 data acquisition device
- USB to Unistat cable
- Dow Syltherm HF heat transfer fluid (~19 L or 5 gallons)
- Safety Siphon
- CryoWedge insulated cover (CryoWedge only)
- Three CryoCassettes (CryoWedge only)
- (10) Celsius®-Paks (S³ Module only)
- (10) Thermocouple holders (S³ Module only)
- Clamp release tool (S³ Module only)

Selecting a Location

Install the system on dry, level surface such as a lab bench or table that supports a total weight of ~400 lbs (~180 kg). The space required for a typical benchtop system is approximately 188 cm (74") wide by 76cm (30") deep. See Figure 5 for proper arrangement and required clearance. The CryoPilot A is air-cooled and requires adequate ventilation. Since the CryoMixer Jr. oscillates back and forth, avoid use of wheeled carts.


Figure 5: Space Requirements

CryoWedge and CryoHose Installation

Caution:

Do not stretch or kink the CryoHoses in order to attach them to the CryoWedge. This may cause leaking. Allow adequate slack.

Caution:

Apply adequate counter torque when tightening the CryoPilot A union fittings in order to avoid damage to the piping inside.

- Place the CryoWedge in the insulated cover provided, making sure that the heat transfer fluid (HTF) inlet and outlet connections are exposed.
- Place the insulated CryoWedge on top of the CryoMixer Jr. Position the CryoWedge such that the motion of the CryoMixer Jr. will allow for mixing along the length of a CryoCassette or from the point of the CryoWedge (equivalent to the position of the core heat exchanger of a CryoVessel) to the outer wall (equivalent to the position of the jacket of a CryoVessel).
- Slide two CryoWraps onto each hose.

Note:

The quick disconnects are factory assembled onto the adapters and hoses. If it becomes necessary to break and then remake connections having tapered threaded fittings, use a thread sealant such as Loctite 567.


Figure 6: CryoHose Installation for CryoWedge

Figure 7: S³ Module Positioning

- Connect the port labeled “Thermofluid out” on the CryoPilot A to the inlet port of the CryoWedge as shown in Figure 6. The inlet port on the CryoWedge is on the outer wall, which is equivalent to the position of the jacket on a CryoVessel.
- Connect the port labeled “Thermofluid in” on the CryoPilot A to the outlet port of the CryoWedge as shown in Figure 6. The outlet port on the CryoWedge is at the point, which is equivalent to the position of the core heat exchanger on a CryoVessel.
- Verify that all hoses are connected properly.

S3 Module and CryoHose Installation

Caution:

Do not stretch or kink the CryoHoses in order to attach them to the S³ Module. This may cause leaking. Allow adequate slack.

Caution:

Apply adequate counter torque when tightening the CryoPilot A union fittings in order to avoid damage to the piping inside.

- Place the S³ Module on top of the CryoMixer Jr. as shown in Figure 7 so that the motion of the CryoMixer Jr. will allow for mixing along the length of a Celsius[®]-Pak.
- Secure the S³ Module to the CryoMixer Jr. plate by tightening the four mounting knobs at the bottom of the module.

- Slide two CryoWraps onto each hose.

Note:

The quick disconnects are factory assembled onto the S³ and hoses. If it becomes necessary to break and then remake connections having tapered threaded fittings, use a thread sealant such as Loctite 567.


Figure 8: CryoHose Installation for S³ Module

- Connect the port labeled “Ther-mofluid out” on the CryoPilot A to the inlet port of the S³ Module. The inlet port is the lower of the two quick-disconnect ports on the S³.
- Connect the port labeled “Ther-mofluid in” on the CryoPilot A to the outlet port of the S³ Module. The outlet port is the upper of the two quick-disconnect ports on the S³.
- Verify that all hoses are connected properly.

CryoPilot A Setup

Warning: All wiring should be performed by a qualified electrician.

- Ensure that the main CryoPilot A power switch is in the “OFF” or “0” position. The switch is red and yellow and is located on the back of the CryoPilot A as shown in Figure 9.

Figure 9: CryoPilot A Main Power Switch

- Have a qualified electrician install a plug (not provided) on the power cord of the CryoPilot A. The plug must match the receptacle and must be rated for the line voltage and current. Refer to System ComponentsSystem ComponentsSystem Components on page 3 for electrical supply specifications.
- Plug the CryoPilot A unit into the electrical supply.
- Turn the main CryoPilot A power switch to the “ON” position.
- Verify that all LEDs light momentarily, and that the rectangular LCD display is illuminated. If not, verify that power is available to the CryoPilot A unit. If power is available to the CryoPilot A unit, but the controller does not come on, please contact Customer Service (see p. 27).

Heat Transfer Fluid Fill

It is recommended that about 19 L (5 gallons) of HTF are on hand for the setup of the benchtop system. Dow Syltherm HF is used as a heat transfer fluid because of its low freezing point and low viscosity over a wide temperature range. For additional information, refer to the Dow Syltherm HF Material Safety Data Sheet.

Note: Dow Syltherm HF is the only approved heat transfer fluid for use in this system.

- Ensure that the drain valve is closed by turning clockwise to the vertical position and closing the lower sealing screw on the side of the CryoPilot A as shown in Figure 10.


Figure 10: CryoPilot A HTF Drain Valve and Sealing Screws

- Remove the fill port cap located on top of the CryoPilot A, behind the sight glass.
- Elevate the container of Syltherm HF above the CryoPilot A.
- Insert the tube end of the Safety Siphon into the CryoPilot A and the valve end into the container of Syltherm HF. Submerge the red valve end in the Syltherm HF and wiggle it up and down to start the siphon. To stop the flow of HTF, lift the red valve end out of the Syltherm HF. Fill the CryoPilot A unit to the 2L mark on the sight glass.
- Start a manual air purge (press Start, Start air purge, 30s interval, OK). Add more HTF as needed so that the sight glass level stabilizes at 2L.
- Verify that there are no leaks in the hose connections.
- Stop the pump (press Stop, Air purge off).
- Secure the CryoWraps over the exposed metal surfaces at each end of the two hoses.

CryoMixer Jr. Setup

- Plug the CryoMixer Jr. into the electrical supply.
- Set the manual operation knob on the front of the CryoMixer Jr. to the "0" position.
- Place the Manual | Automatic switch on the back of the CryoMixer Jr. to the manual position (down) as shown in Figure 11.

Figure 11: Rear of CryoMixer Jr.

- Turn the speed knob to the fully counterclockwise (0 RPM) position.
- Turn the manual operation knob to the "1" position.
- Slowly turn the speed knob in the clockwise direction until the proper RPM is achieved and record displayed RPM as read on the LCD display on the CryoMixer Jr. The standard settings are listed below:
 - CW20 and CW125 CryoWedges: 45 RPM
 - CW300 (30") and CW300 (34") CryoWedges: 35 RPM
 - S³ Module: 120 RPM

Note:

Mixer speed settings may vary according to customer product characterization. Contact Sartorius Stedim Biotech to discuss variations from the above recommendations.

- Place the Manual | Automatic switch on the back of the CryoMixer Jr. to the automatic position (up).

Computer Setup

- Plug the computer and display into the electrical supply. It is recommended that you use a surge protector (not included).
- Turn on the computer and display. Answer the questions in the Windows "first boot" setup dialogs. It may be necessary to obtain some information from your Information Technology (IT) department.

Caution:

Computers supplied by Sartorius Stedim Biotech do not have antivirus software installed. It is recommended that you install antivirus software according to your company's policy.

Data Acquisition Setup

- Connect the National Instruments USB-9213 to the computer with a USB cable.
- Connect one end of the USB to CryoMixer Jr. cable to the receptacle at the rear of the CryoMixer Jr. Plug the other end into a USB port on the computer.

CryoPilot Software

- For S³ Modules, plug the bundle of male thermocouple connectors extending from the S³ Module into the thermocouple breakout cable. Plug connector 1 into receptacle 1, connector 2 into receptacle 2, and so on. Plug the Celsius[®]-Pak thermocouples into the receptacle strip located within the S³ hood enclosure area.
- Connect one end of the USB to Unistat cable to a USB port on the computer and the other end to the serial port on the underside of the ComBox as shown in Figure 12. The ComBox is located at the rear of the CryoPilot A.

Figure 12: CryoPilot A Serial Port

Overview

The benchtop system is controlled and monitored by CryoPilot software, which is a graphical application for computers with Microsoft Windows. The CryoPilot software provides users with flexible execution of freeze and thaw operations.

The following diagram illustrates the relationships between system components. The CryoPilot software interfaces with the CryoPilot A chiller and the CryoMixer Jr. shaker via USB, and with up to 16 thermocouples through a National Instruments data acquisition device.

Figure 13: CryoPilot System

System Requirements

CryoPilot software is factory-installed on computers ordered from Sartorius Stedim Biotech as part of a benchtop system. Alternately it is available on CD-ROM for systems that meet the following requirements:

- Operating System: Microsoft Windows XP Professional, 32-bit or Microsoft Windows 7 Professional, 32-bit.
- Data Acquisition Device: National Instruments USB-9213. National Instruments PCI-4351, USB-4350 (discontinued), and DAQCard-4350 (discontinued) are only supported on Windows XP.

CryoPilot software is packaged with an installer that will automatically install the CryoPilot software and all required third-party drivers on your computer. You must have administrative privileges to install CryoPilot.

Components

By default, CryoPilot software is installed in C:\Program Files\Sartorius Stedim Biotech\CryoPilot 5. This will hereafter be called the installation directory. The following is a description of the files and folders contained in the installation directory:

Component	Description
Profiles/	Directory for standard profile files supplied with CryoPilot .
BouncyCastle.Crypto.dll	Support library
CryoPilot.chm	Help file
CryoPilot.exe	The executable for CryoPilot software. Shortcuts to this file are placed on the desktop and in the Start menu.
CryoPilot.exe.config	Settings file
ICSharpCode.SharpZipLib.dll	Support library
libftdi.dll	Support library
log4net.dll	Support library
NationalInstruments.Common.dll	Support library
NationalInstruments.Common.Native.dll	Support library
NationalInstruments.DAQmx.dll	Support library
NPlot.dll	Support library
unins000.dat	Uninstaller support file
unins000.exe	Uninstaller

In addition to the files in the installation directory, the CryoPilot installer creates a folder named CryoPilot in each user's My Documents folder. This will hereafter be called the user directory. The following is a description of the files and folders contained in the user directory:

Component	Description
Data/	Default directory for data files saved by users. Users may opt to store files elsewhere on the computer or on a network drive.
Custom Profiles/	Default directory for custom profiles saved by users. Users may opt to store files elsewhere on the computer or on a network drive.

Files that are subject to change but common to all users are stored in C:\Documents and Settings\All Users\Application Data\Sartorius Stedim Biotech\CryoPilot 5 on Windows XP and in C:\ProgramData\Sartorius Stedim Biotech\CryoPilot 5 on Windows 7.

Component	Description
log.txt	File used for troubleshooting software
settings.xml	Settings file

The CryoPilot installer also installs several additional libraries: Microsoft .NET Framework 2.0, National Instruments NI-DAQmx, and USB drivers. The installer also provides optional installation of National Instruments Traditional NI-DAQ (Legacy) for support of legacy data acquisition devices including PCI-4351, USB-4350, and DAQCard-4350.

Installation

To install CryoPilot software, insert the CD-ROM and double-click on setup.exe. Setup requires administrative privileges; contact your system administrator if you do not have sufficient privileges.

The default settings are generally adequate. If you are using legacy hardware such as National Instruments PCI-4351, USB-4350, or DAQCard-4350, be sure to check the box "Support for PCMCIA-4350, USB-4350, PCI-4351".

During the setup process you will be prompted twice about installation of unsigned drivers. Be sure to accept the dialogs and continue.

Once the setup process is complete, you will be prompted to restart your computer.

It is critical that the computer be configured so that it does not suspend or power down while a freeze/thaw process is running.

For Windows 7, go to Control Panel > System and Security > Power Options. Select High performance profile. Click Change plan settings, then Change advanced power settings. Set Hard disk > Turn off hard disk after to 0 (never). Set USB settings > USB selective suspend setting to Disabled. Click OK, then close the Control Panel window.

For Windows XP, go to the Control Panel > Power Options. Select the Presentation power scheme. Click OK, then close the Control Panel window.

User Interface

Main Window

The main window is the primary screen used in the CryoPilot software. This window contains menu items to interact with the program, a chart to graphically display process data, and a legend to select the traces to display on the chart.


Figure 14: Main Window

ID	Control	Description
1	Title bar	Displays the name of the application and the name of the current data file, if any
2	Menu bar	Contains menus for performing common tasks
3	Toolbar	Contains buttons for performing common tasks
4	Status bar	Shows the status of the equipment as Online or Offline
5	Chart	Displays a plot of temperature vs. elapsed time
6	Elapsed time indicator	Shows how much time has elapsed since the process started
7	HTF Setpoint legend item	Displays the heat transfer fluid (HTF) setpoint
8	HTF Supply legend item	Displays the actual HTF temperature
9	Product N legend item	Displays the temperature of the Nth type T thermocouple. Up to 14 thermocouples may be monitored, depending on hardware.

Start Process Dialog

The Start Process dialog prompts the user for information that will be used during the process.

Figure 15: Start Process Dialog

ID	Control	Description
1	Chart title	The title shown at the top of the chart. The default is "Temperature vs. Elapsed Time".
2	Data filename	The complete file name where data for this process will be saved. The default is [User Directory]\Data\[date code].cpd, where [date code] is a numeric string composed of year, month, day, hour, minute, second.
3	Data filename selector	Opens a dialog to select a folder and filename where the data for this process will be saved.
4	Profile type selector	Gives three options: None, for no profile; Standard, for a Sartorius Stedim Biotech – supplied standard profile, or Custom, for a user-supplied profile.
5	Application selector	Gives four options when Standard profile type is selected: – S ³ 30mL – for S ³ systems with ten 30mL Celsius®-Paks – S ³ 100mL – for S ³ systems with ten 100mL Celsius®-Paks – CryoWedge 30" – CryoWedge 34"
6	Process selector	Gives three options when Standard profile type is selected: – Freeze – performs freeze only – Thaw – performs thaw only – Freeze and Thaw – performs freeze followed by thaw
7	Profile filename	The complete file name where the custom profile can be found
8	Profile filename selector	Opens a dialog to select a folder and filename where the custom profile can be found.
9	Description	A field for entry of additional information about the process.
10	OK button	Accepts the information, closes the dialog, and starts the process.
11	Cancel button	Ignores the information and closes the dialog without starting the process.

Options Dialog

CryoPilot software may be configured using the CryoPilot Configuration dialog. From the main window, select the Tools > Options menu item.


Figure 16: Options Dialog – Certificate Tab

ID	Control	Description
1	Certificates list	Displays a list of available digital certificates. For each certificate the list shows who issued the certificate, to whom the certificate was issued, the friendly name, and the expiration date.
2	Create button	Creates a new self-signed certificate and installs it in the Trusted Root Certification Authorities folder.
3	View button	Views additional details for the selected certificate.
4	Manage button	Opens the Microsoft Management Console for additional certificate management tasks such as certificate renewal and deletion.
5	OK button	Accepts the changes and closes the dialog.
6	Cancel button	Discards the changes and closes the dialog.
7	Hardware tab	Switches to the hardware tab.

Figure 17: Options Dialog – Hardware Tab

ID	Control	Description
1	Data acquisition device selector	The “None” option configures the application for no data acquisition. If a compatible National Instruments data acquisition device is detected, it will be automatically selected.
2	Chiller selector	Selecting the “None” option configures the application for no chiller. Selecting the “Huber Unistat” option configures the application to collect data from a Huber chiller equipped with a Unistat controller, such as a CryoPilot A.
3	COM port selector	Shows a list of available serial (COM) ports.
4	Shaker selector	The “None” option configures the application for no shaker. If a USB to CryoMixer Jr. cable is detected, it will be automatically selected.
5	OK button	Accepts the changes and closes the dialog.
6	Cancel button	Discards the changes and closes the dialog.
7	Certificates tab	Switches to the Certificates tab.

Getting Started

Before a profile can be started the software must be configured. Make sure that the data acquisition device (DAQ) and chiller are connected and powered up. Then select the Tools > Options menu item to display the Options dialog.

Selecting a Digital Certificate

First, a digital certificate must be selected. CryoPilot uses a digital certificate to guarantee the authenticity and integrity of saved data files. The Certificates tab shows a list of available certificates. It is ideal to use a certificate issued by a trusted certification authority (CA). If your information technology (IT) organization already has a digital certificate infrastructure established then they may be able to provide a certificate for use with CryoPilot.

In the event that you do not have a certificate to use with CryoPilot, you can create what is known as a self-signed certificate. "Self-signed" means that the person who creates the certificate also vouches for its legitimacy. For Click the Create button, and then click the Yes button on the dialog that appears. The newly created certificate will be shown in the list of available certificates.

Configuring Hardware

Next, the hardware must be selected. The Hardware tab contains settings for common hardware configurations. The most common arrangement consists of a data acquisition device (DAQ), a chiller, and a shaker. In most cases the hardware will be detected automatically.

If you will not be running any processes and will only be using CryoPilot to review data files, then select "None" for all Data Acquisition, Chiller, and Shaker.

To save the options, click OK.

Program Flow

When CryoPilot software is launched, a splash screen appears briefly and then the main window is displayed.

At this point the software is idle and will not do anything until the user starts a process.

To start a process, select the Process > Start menu item or click the Start Process button on the toolbar. The Start Process dialog will be displayed. Once the user clicks OK, the software attempts to connect to the data acquisition device, the shaker, and the chiller, if configured. If any of these devices is configured but not connected then the user is prompted to check connections and retry or cancel.

Once a process is started, the temperature chart and the legend are reconfigured to display the appropriate plots, communications are established with the chiller, and the data acquisition device (DAQ) is initialized. Once initialization is complete, the initialization message disappears and first set of data points is plotted. Subsequent data points are displayed at the sampling interval, which is fixed at 15 seconds.

When configured to use a chiller, the CryoPilot software uses a serial link to send commands to and receive responses from the chiller.

During a process, the connection status is displayed in the lower left corner of the main window. A green "Online" message indicates a good connection with the data acquisition device and the chiller. A red "Offline" message indicates that something is wrong with the connection.

The Stop Process button ends data acquisition and turns off the chiller and mixer. The user may then review the data. The data file must be closed before the user can start another process.

Reviewing Data

CryoPilot can open and display files saved in CryoPilot 4 and CryoPilot 5. Note that files created by CryoPilot 4 are supported for review purposes only. All new data files are saved in CryoPilot 5 format so that data integrity and authenticity can be verified via digital signature.

Only one data file can be open at a time.

To open a data file, do one of the following:

- Click the Open button on the toolbar
- Open the File menu and click Open
- Use the keyboard shortcut Ctrl+O

Using the Open Data File dialog, select the file you wish to open. Be sure to select the type of file in the "Files of type" dropdown.

CryoPilot 5 uses digital signatures to prove that the data has not been modified since it was originally saved. If the integrity and authenticity of the data file cannot be verified by means of a digital signature then you will be notified and prompted to continue or cancel. There are two common cases where this notification will appear:

- When opening a CryoPilot 4 data file. CryoPilot 4 files do not contain a digital signature.
- If the certificate used to sign the data file is not trusted.

Once the data file is loaded, the chart will be populated with the data from the data file, and the data on the chart can be explored just as if the data was coming from a live process.

By right-clicking on the chart and clicking on the Review item, you can scroll right and left on the chart and verify the values for each trace at a given point in time. As the review line is scrolled right or left, the legend is updated to show the timestamp and trace values for that point in time.

Using the Chart

The chart is a graphical display of temperature vs. elapsed time for a run. Elapsed time is displayed on the X axis in units of hours and temperature is displayed on the Y axis in units of degrees C. A legend is provided to identify and show the current value for each plot.

Panning and Zooming

By default, chart scales are automatically adjusted to fit the available data. This auto-scaling feature may be turned on or off using the Auto Sscale item X and Autoscale Y options on the pop-up menu obtained by right-clicking on the chart.


Figure 18: Chart Pop-up Menu

The scale extents can be changed directly by dragging the X and Y axes. The chart display can be panned by dragging the chart with the left mouse button held down. When changing scale or panning the chart, Auto Scale will be turned off. Double-click the chart to turn Auto Scale back on.

Legend

The legend identifies each of the plots on the chart. Each plot has a legend item showing its name, color, units, and current value. To hide a plot, double-click its legend item; it will appear dimmed. To show the plot, double-click the legend item again. This can be useful when trying to look at data files with many plots lying on top of each other.

Review Mode

Review mode allows the user to review data at any point during the process. Review mode is turned on and off using the pop-up menu shown in Figure 18. When in review mode, a vertical line is displayed on the chart, and the elapsed time and legend items display values that correspond to that location on the chart.

Exporting Data

CryoPilot provides two different methods to export data. It can export an image of the current chart view or it can export raw data for viewing in a spreadsheet application.

Export Chart

To export the chart, select the Edit > Copy > Chart menu item. The chart is copied to the Windows clipboard as a bitmap image. The chart may be pasted to other applications by selecting Edit>Paste or pressing Ctrl+V. This option is useful where a quick snapshot of the data is needed, for example in a slide presentation. Use the Copy>Data option if the data must be further analyzed.

Export Data

To export the data, select the Edit > Copy > Data menu item. The raw data is copied to the Windows clipboard in tab-delimited format. The data may be pasted to spreadsheet applications such as Microsoft Excel or OpenOffice Calc by selecting the destination cell and selecting Edit>Paste or pressing Ctrl+V.

Printing

CryoPilot can create a report for the open data file. The report is a snapshot of the current chart view including the title and legend. The report may be previewed on-screen or sent to a printer. If you have a PDF print driver such as Adobe Acrobat installed on your computer, you can save the report to a PDF file on disk.

To print a report, do one of the following:

- Click the Print button on the toolbar
- Open the File menu and click Print
- Use the keyboard shortcut Ctrl+P.

The Print dialog will appear, prompting you to select the destination printer and the number of copies. Other settings such as margins and paper size may be set in Page Setup. Click OK to send the report to the selected printer.

Profiles

CryoPilot software controls the chiller and the mixer by following a profile. A profile is a recipe with instructions to change the chiller temperature setpoint and the mixer mode setpoint (on|off) as a function of time. CryoPilot software is supplied with a set of standard profiles that are suitable for many S³ Module and CryoWedge applications.

It may be necessary to create a new profile to suit a specific application. When designing a new profile, there are several factors to consider:

- Which operations are to be performed – freeze, thaw, or both?
- Are multiple freeze and thaw cycles required?
- What is the initial temperature of the material to be frozen or thawed? If this is variable, an equilibration step at the beginning of the profile can improve reproducibility.
- Should the equipment remain running or turn off after the freeze or thaw is complete?

Adjustment of the freeze and thaw rates is performed empirically. If a lower endpoint is desired for a freeze, it may be necessary to increase the duration of the final setpoint or reduce the the temperature adjusted. When trying to thaw quickly without overshoot, use a high temperature setpoint for a short duration and then lower the temperature setpoint before the target temperature is exceeded.

Profile Editing

Profiles may be created and edited manually with text editors (e.g. Microsoft Notepad) or spreadsheet programs (e.g. Microsoft Excel). When using a text editor, always use the tab key to separate columns; do not use the space key. For a new line, use the Return or Enter key after entering the status for the Mixer.


Figure 19: A profile in Microsoft Notepad.

If you use a standard profile as a starting point, be sure to save it with a new file name and then make changes.

Last Step Behavior

By default, once the CryoPilot software executes the last step of a profile it will leave the chiller and mixer as-is until the user stops the profile by pressing the STOP button. In some cases it is useful to turn the chiller and mixer off at a certain elapsed time. This is done by placing a TURN CHILLER OFF instruction after the STOP instruction. Note that this cannot be done using the CryoPilot profile editor; it must be done manually with a text editor or spreadsheet application. For example, the following profile will turn the chiller and mixer off at 311 minutes:

```
0      -40    OFF
310.00 -60    OFF
STOP
311.00  TURN CHILLER OFF
```

The TURN CHILLER OFF instruction will be visible in the CryoPilot profile editor, but will have no temperature associated with it.

File Format

Profiles are stored on disk as tab-delimited data files. The character encoding is ASCII and the file format is case-insensitive. All lines found before the START token are a description of the file. This means that the description cannot contain the word "start". Between the START token and the optional STOP token, each line represents a step in the profile. The start time identifies when the step should start; it has units of elapsed minutes and is expressed in decimal notation. The start time of the first step must be 0, and start times for successive steps must increase monotonically. Note that the start time for each step references the profile start time, so three 30-minute steps would have start times of 0, 30, and 60. The chiller setpoint has units of °C and is expressed in decimal notation. If a STOP token is present, the next line may be a TURN CHILLER OFF token.

The following syntax is used to describe the file format:

```
\n      Newline character, here defined as CRLF
\s      Zero or more whitespace characters
\t      Tab character
<variable> Variable element
[FOO]   Optional element
OFF|ON  Enumerated value, element must be one or the other
...     Indicates repeating structure
```

Format:

```
[<description>\n]
START\s\n
<start time>\t<chiller setpoint>\t<mixer = OFF|ON>\n
.
.
.
[STOP\n
 [<start time>\tTURN CHILLER OFF]]
```

Example (implied \n):

This is a test profile with three steps.

```
START
0      -20    OFF
30     -60    ON
200    10     OFF
```

Data Files

During a run, the CryoPilot software monitors the benchtop system components and records key process values in a data file on the computer.

CryoPilot 5 software can open data files created by CryoPilot 4 as well as data files created by CryoPilot 5. However, all new data files saved by CryoPilot 5 will be saved in CryoPilot 5 format. This data file contains the following process values:

- Chiller heat transfer fluid temperature setpoint and actual heat transfer fluid temperature
- Chiller mode setpoint (On | Off)
- Mixer mode setpoint (On | Off)
- Temperature measured by up to 14 Type T thermocouples

By default, CryoPilot 5 data files have a .cpd file extension.

A CryoPilot data file is a ZIP file containing two files:

- data.txt, which is a tab-delimited data file
- data.txt.p7s, which is a digital signature in PKCS#7 format.

Spreadsheet programs such as Microsoft Excel can import the tab-delimited data file, but it is easier to transfer the data using the File > Copy > Data menu item. Do not attempt to open the data file for a run in progress as the spreadsheet application may lock the file, thus preventing the CryoPilot software from writing additional data.

File Format

The data.txt file and the results of the File > Copy > Data command both contain tab-delimited text. The character encoding is ASCII and the file format is case-insensitive. The first line must contain the creator and version as specified below. All temperature values have units of °C.

The following syntax is used to describe the file format:

\n Newline character, here defined as CRLF
\s Zero or more whitespace characters
\t Tab character
<variable> Variable element
[FOO] Optional element
OFF|ON Enumerated value, element must be one or the other
... Indicates repeating structure

Note that times are specified as Universal Coordinated Time (UCT), which may differ from local time depending on your time zone.

Where no data is available, the file will contain the symbol "NaN" (Not a Number).

Format:

```
CryoPilot 5.0\n
METADATA\n
Title\t<title>\n
Process Start\t<timestamp format=yyyy-MM-dd HH:mm:ssZ>\n
Computer\t<computer name>\n
User\t<user name>\n
Description\t<description>\n
START\n
Date/Time(UTC)\tElapsed Time(h)\tChiller Status\tHTF Supply Temperature Setpoint(C)\tHTF Supply Temperature(C)\tMixer Sta-
tus\tProduct Temperature 1(C)\tProduct Temperature 2(C)\t...Product Temperature n(C)\n
<timestamp format=yyyy-MM-dd HH:mm:ssZ>\t<elapsedhours format=0.000000>\t<chiller = Off|On|Unknown>\t<value for-
mat=0.0>\t<value format=0.0>\t<mixer = Off|On|Unknown>\t<value format= 0.0>\t<value format= 0.0>\t...<value format=0.0>\n
.
.
.
.
```

Example with four thermocouples (implied \n):

```
CryoPilot 5.0
METADATA
Title           Temperature vs. Elapsed Time
Process Start   2009-03-03 21:49:22Z
Computer        DEMO_LAPTOP
User            DEMO_LAPTOP\Administrator
Description
START
Date/Time(UTC)  Elapsed Time(h)  Chiller Status  HTF Supply Temperature
Setpoint(C)     HTF Supply Temperature(C)  Mixer Status    Product 1
Temperature(C)  Product 2 Temperature(C)    Product 3 Temperature(C)
Product 4 Temperature(C)
2009-03-03 21:49:31Z  0.002500  On  0.0  21.0  On  22.8  25.4  NaN  NaN
2009-03-03 21:49:46Z  0.006667  On  0.0  20.7  On  22.8  25.3  NaN  NaN
2009-03-03 21:50:01Z  0.010833  On  0.0  20.3  On  22.8  25.3  NaN  NaN
2009-03-03 21:50:16Z  0.015000  On  0.0  20.0  On  22.8  25.2  NaN  NaN
```

Diagnostics | Troubleshooting

Error and status messages are displayed in the message area at the lower right hand corner of the main window. This section describes common errors and what may be done in response.

Chiller

In the event of problems with the chiller, the following steps may be helpful:

- Check both ends of the cable between the computer and the chiller.
- Check the chiller's error log. For more information, refer to the Huber Unistat 705 Operating Instructions.
- Check that the chiller communication settings are configured properly. The CryoPilot A should be set to use baud rate 9600 and protocol RS-232. For more information, refer to the Huber Unistat 705 Operation Instructions.
- As a last resort, exit the CryoPilot software, cycle power to the chiller, and restart the CryoPilot software. Try to ascertain the cause of the problem before starting another run. Make sure that the settings in the Options dialog are correct.

If a power failure occurs during a run the chiller may get disconnected or the data from the chiller may be corrupted. The presence of NaN (not a number) in the temperature fields of the data file is another indication of a chiller error.

Shaker

If the shaker does not move when it should, the following may be helpful:

- Check the cable between the computer and the shaker.
- Verify that the mode switch at the rear of the mixer is in the Auto position.
- Verify that the speed setting is correct.

Data Acquisition Device | Thermocouples

The following are common errors related to the data acquisition hardware or the thermocouples:

- A sensor does not respond to changes in temperature as the other sensors do.

Should one of the above occur, the following hints may be helpful:

- Verify that the sensor in question is plugged in to the proper receptacle in the thermocouple breakout cable attached to the data acquisition device.
- Be sure that each sensor is properly located and secured. A small error in location can cause large deviations compared to other sensors.
- Check connections between the computer and the data acquisition device.
- A sensor may be damaged. Check the sensor with a calibrated reader (e.g. Fluke 725) capable of measuring Type T thermocouples.

21 CFR Part 11

CryoPilot 5 software has been designed to facilitate compliance with the requirements of 21 CFR Part 11.

Digital signatures are used to prove that the authenticity and integrity of CryoPilot 5 data files, even on open (i.e. untrusted) systems. The terms “digital signature” and “electronic signature” should not be confused. A digital signature is a cryptographic method, whereas an electronic signature is an electronic equivalent of a handwritten signature with legal implications. References to electronic signatures in 21 CFR Part 11 do not apply to the digital signatures used in CryoPilot 5 data files.

Some general notes about implementation:

- CryoPilot 5 relies on the operating system for account management, authentication, and authorization functions.
- The digital signature algorithm used by CryoPilot 5 is SHA-1 with RSA; it is written in PKCS#7 format.
- The file system is treated as an untrusted system.
- When a process is stopped, the data file is closed and cannot be opened again for modification. Any alteration of the data by other means, whether accidental or malicious, will be detected by the digital signature mechanism.
- Once data is exported from CryoPilot (e.g. to a spreadsheet), its authenticity and integrity can no longer be ensured.

System Operation

Preparing for a Run

CryoWedge

To prepare the CryoWedge for use:

- Remove the upper half of the CryoWedge insulated cover.
- Remove the wing nuts holding the transparent cover and set aside.
- Remove the transparent cover. Fill the CryoWedge with the test solution using a graduated cylinder. Standard volumes for each size CryoWedge are as follows:

CryoWedge Size	Volume (mL)
CW20 (12 inch)	350
CW125 (20 inch)	650
CW300 (30 inch)	3000
CW300 (34 inch)	4000

- Replace the transparent cover.
- Replace the wing nuts and tighten them hand tight. Do not over-torque or strip the threads.
- Replace the upper half of the CryoWedge insulated cover.
- Ensure that the CryoMixer Jr. cable is plugged into the CryoMixer Jr. and the computer.
- Ensure that at least six thermocouples are plugged into the data acquisition device, starting with position 1.
- Place five of the thermocouples through the pre-drilled holes in the transparent cover. Secure the thermocouples 2" below the surface of the test solution using plastic spacers or cut pieces of silicon tubing.
- Place the remaining thermocouple on the top of the transparent cover as a reference temperature at the outer surface of the CryoWedge.
- Record the number and position of each thermocouple for later reference.

S³ Module

Celsius[®]-Paks are delivered pre-sterilized by gamma irradiation. They are supplied with female Luer fittings and caps on 2 fill | drain ports. A variant with C-Flex tubing extensions is also available for sterile filling. All Celsius[®]-Paks have a closed silicone/EVA thermowell for insertion of a thermocouple probe.

In order to obtain consistent results, always load the S³ Module with the same number of Celsius[®]-Paks. Ensure that each Celsius[®]-Pak is filled to the same volume. The thermocouple should be placed at the same depth every time.

The standard freeze profiles provided with the CryoPilot software are designed for 10 Celsius[®]-Paks, either 30 or 100 mL, filled to nominal volume. Place thermocouples at a minimum depth of 1 cm below the liquid surface of Celsius[®]-Paks. Blank Celsius[®]-Paks should be filled with sample buffer. Do not use deionized water in blanks as it typically supercools. If you use fewer than 10 Celsius[®]-Paks, the freeze completes more rapidly and the profile does not transfer directly to a larger scale.

Step by step procedure for Celsius[®]-Pak prep, filling and loading:

- Prepare sample and | or buffer for filling Celsius[®]-Paks.
- Insert label in label holder at the bottom of Celsius[®]-Pak. A label with dimensions of ~2 × 5 cm fits entirely within the label holder. It is recommended to use waterproof material and ink for the label. Do not write on the Celsius[®]-Paks as the solvents used in the ink may penetrate the Celsius[®]-Pak and contaminate the product.

Use a syringe or peristaltic pump to fill Celsius[®]-Paks with sample or buffer through one of the fill | empty ports. Connect the syringe or pump via a male Luer[®] fitting or welded with the optional C-Flex line. Fill Celsius[®]-Paks with a maximum volume of 30 mL for the 30 mL Celsius[®]-Pak or 100 mL for the 100 mL Celsius[®]-Pak.

Caution:

Do not overfill Celsius[®]-Paks or purge all air from the top of the Celsius[®]-Paks as this will cause sample and ice expansion into port tubing.

Close the fill tube with the supplied Luer[®] cap or pinch clamp. Wipe any excess moisture from all surfaces inside the S³ Module, prior to starting a freeze.

Celsius[®]-Paks are held in the S³ Module by a spring clamp mechanism that holds the Celsius[®]-Pak and acts as an active heat transfer surface. Perform the following steps for loading the S³ Module:

Attach thermocouple holder to Celsius[®]-Pak by aligning the female Luer lock on thermocouple holder with the male Luer lock on the top of the thermowell. After fastening the Luer lock to the thermowell, note that the front of the thermocouple holder should face the same direction as the front label of the Celsius[®]-Pak.

Figure 20:
Thermocouple Holder Attached to Celsius[®]-Pak

- Measure 1 cm from the end of the thermocouple and place a mark at this point for placement of thermocouple within the Celsius[®]-Pak.
- Slowly unscrew the thermocouple locking collar, taking care to collect the two small o-rings within the locking collar. Place the locking collar and two o-rings over the end of the thermocouple and screw down the locking nut just far enough to ensure that it will not fall off the thermocouple holder.
- Lower the thermocouple into thermowell. Never push hard enough to pierce through the end of the thermowell. For accurate thermal monitoring, always place the thermocouples 1 cm below the surface of the liquid in the Celsius[®]-Pak by ensuring that the mark on the thermocouple meets the liquid surface.

Figure 21: Thermocouple 1 cm below Liquid Level Surface

Tighten the locking collar to secure the thermocouple in place within the thermowell.

Stand facing the left side of the S³ Module; insert a Celsius[®]-Pak between the clamp mechanisms, starting at slot 10; the slot furthest from the HTF supply and return. Open the mechanism by placing the clamp release tool over the 2 guide pegs on the 2 stainless steel pinch blocks for the slot and pressing down.

Figure 22: Opening the Pinch Clamp Mechanism

- While pressing the clamp release tool down, Slide each of the sides of the Celsius[®]-Pak down between a black block and a stainless steel pinch block of the clamping mechanism. Position the top of the Celsius[®]-Pak so it is even with the top of the stationary black metal bar of the clamp mechanism. Note that the Celsius[®]-Pak must face left in order to fit properly into the S³ Module.

Figure 24: Celsius[®]-Pak Placement and Thermocouple Connection

- Repeat this procedure from slot 9 to 1. Figure contains a diagram of the final placement of all Celsius[®]-Paks and thermocouples.

Note:

For reproducibility, load the S³ Module with the same number of Celsius[®]-Paks, fill each Celsius[®]-Pak with the same amount of sample or buffer and place the thermocouples in the thermowells 1 cm below the surface of the liquid. All freeze profiles developed at Sartorius Stedim Biotech perform correctly only with a fully loaded system (10 Celsius[®]-Paks). For that reason, results achieved with the freeze profiles provided with the installation software are directly scalable to 16.6 Liters. Changing the number of Celsius[®]-Paks affects the scalability of these profiles.

Figure 23: Thermocouple Holder and Celsius[®]-Pak Installation

- Withdraw the clamp release tool. The Celsius[®]-Pak will be held tightly on each side between two metal blocks, one black and one stainless steel. Connect the male thermocouple plugs to the corresponding thermocouple receptacle strip.

- Diagram the position and number of each thermocouple in a lab notebook or in accordance with company policy.
- Close and lock the cover of the S³ Module.


Figure 25: Thermocouple and Celsius[®]-Pak placement in S³ Module.

Running a Profile

Note:

When thawing with the S³ Module, ensure that the 4 S³ Module mounting knobs are tight and that the module is held securely on the CryoMixer Jr. plate.

Caution:

Never exceed a value of 120 RPM for the CryoMixer Jr.

- Double click the CryoPilot software icon to initialize the CryoPilot Software.
- Select the Process > Start menu item to display the Start Process dialog.
- Provide a file name for the data file where time and temperature data will be stored for the new process.
- Select the Standard profile type and select your application and desired operation (freeze, thaw, or freeze and thaw). Alternately, select a custom profile.
- Click the OK button to start the process.
- Verify that the chiller responds to the commands from the software. The status indicator in the lower left hand corner should read "Online". There will be a 2-5 minute delay before the compressor starts and temperature regulation commences.
- Allow the profile to run to the end of the programmed time.
- Select the Process > Stop menu item to stop the process.

Removing Frozen Celsius[®]-Paks from the S3 Module

In some cases it is desirable to remove Celsius[®]-Paks from the S³ Module while they are still frozen.

Remove the transparent top cover of the S³ Module.

Note:

If the unit is left open for an extended period of time while maintained at low temperature, frost builds up on the clamps, Celsius[®]-Paks and thermocouple holders making Celsius[®]-Pak removal difficult. Once the unit is open, quickly remove the frozen Celsius[®]-Paks and then replace the hood.

- Remove the thermocouples from the Celsius[®]-Paks. If a thermocouple does not readily pull out, gently twist until it pulls free. Remove the thermocouple holder from each slot. If condensation builds up on the surfaces of the clamps, it forms an 'ice bond' between the thermocouple holders and the black metal cooling blocks. This makes removal of the thermocouple holders nearly impossible. Wipe excess moisture from all surfaces prior to starting a freeze.

- Remove the frozen Celsius[®]-Paks from the S³ Module using the clamp release tool to slide the Celsius[®]-Paks up and out. Place in an insulated transport container.

- In the freezer, store frozen Celsius[®]-Paks in a protective container at the desired temperature. If performing stability tests or formulation studies, it is best to use a controlled freezer with thermal monitoring.

Dispensing Thawed Samples from the S3 Module

- Remove thermocouples and thermocouples holders.
- Remove the thawed Celsius[®]-Paks from the S³ Module.
- Invert each Celsius[®]-Pak three times to ensure uniform product mixture.
- Remove desired amount of sample from Celsius[®]-Pak via the unused fill | empty port. This can be done by syringe, pump, or gravity drain.
- Discard used Celsius[®]-Paks. Since Celsius[®]-Paks are sometimes deformed during freezing and thawing, upon reuse they do not fit as well between the clamping mechanisms. This can cause prolonged freezing and thawing due to poor heat transfer.

System Maintenance

Maintenance may be performed by the end-user or by Sartorius Stedim Biotech. Contact the Sartorius Stedim Biotech Service Department to establish a Preventive Maintenance Contract.

In the event of a problem, first check the error code section of the Huber Unistat 705 Operating Instructions. If you are unable to remedy the situation, call Customer Service (see p. 27).

Cleaning

Exposed equipment surfaces are generally compatible with spray and wipe with isopropyl alcohol. Stainless steel surfaces may optionally be cleaned with stainless steel polish.

The internal cabinet of the S³ contains many parts made of anodized aluminum which is incompatible with basic solutions such as sodium hydroxide (NaOH) or oxidizers. These materials should not be used for cleaning the internal cabinet. Rather, isopropyl alcohol or other non oxidizing cleaning solutions are recommended if cleaning is necessary.

Changing the Heat Transfer Fluid

Change the heat transfer fluid if it becomes visibly contaminated with water or particulate, or if you notice a decline in heat transfer efficiency (i.e. when the unit does not reach or maintain set point temperature). Always change out the entire amount of HTF. If there is no improvement in heat transfer efficiency after changing the HTF, call Customer Service (see p. 27).

To change the HTF, use the following procedure:

- Verify that the system HTF drain valve located on the front of the unit is closed. The slot in the screw cap is vertically aligned when closed. See Figure 10: CryoPilot A HTF Drain Valve and Sealing Screws, p. 10.
- Remove the sealing screw in the protruding stainless tube on the left side of the unit and connect a piece of flexible tubing over the end of tube. Attach securely with a hose clamp.
- Put the far end of the tubing in a container that will hold at least 12 L (3 gallons) of HTF.
- Open the HTF drain valve by turning counter-clockwise until the slot is horizontally aligned.
- If possible, attempt to lift CryoWedge or S² module to drain residual HTF from the application module back down into the CryoPilot A.
- When HTF stops flowing, close the HTF drain valve, remove the tubing from the sealing screw tube, and replace the sealing screw.
- Fill the CryoPilot A system with new HTF according to the procedure in Heat Transfer Fluid Fill on p. 9.

Information and Instructions on Disposal and Repairs

Packaging that is no longer required must be disposed of at the local waste disposal facility. The packaging is made of environmentally friendly

materials that can be used as secondary raw materials.

The equipment, including accessories and batteries, does not belong in your regular household waste. EU legislation requires its Member States to collect electrical and electronic equipment and dispose of it separately from other unsorted municipal waste with the aim of recycling it.

In Germany and many other countries, Sartorius AG takes care of the return and legally compliant disposal of its electrical and electronic equipment on its own. These products may not be placed with the household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other Member States of the European Economic Area (EEA), please contact our service technicians on location or our Service Center in Goettingen, Germany:

Sartorius AG
Service Center
Weender Landstrasse 94-108
37075 Goettingen, Germany

In countries that are not members of the European Economic Area (EEA) or where no Sartorius affiliates, subsidiaries, dealer, or distributors are located, please contact your local authorities or a commercial disposal operator.

Prior to disposal and/or scrapping of the equipment, any batteries should be removed and disposed of in local collection boxes.

Sartorius AG, its affiliates, subsidiaries, dealers, and distributors will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the accompanying manual or visit our Internet website (www.sartorius.com) for comprehensive information that includes our service addresses to contact if you plan to send your equipment in for repairs or proper disposal.

Contact Information

24 Hour Customer Service
+1.866.424.5600
(toll free North America only)
+1.707.747.5614
(outside North America)

EC Declaration of Conformity

Company	Sartorius Stedim Systems GmbH	
Address	Schwalzenberger Weg 23-29 34212 Mitteln, Germany Phone +49 551 30833, Fax +49 551 3083769 www.sartorius-stedim.com	
	<p>We hereby declare that based on the design, construction and version placed on the market, the device designated below fulfills the relevant fundamental safety requirements and health regulations specified by the pertinent EC Directive.</p> <p>This Declaration shall become invalid should any modifications be made to the device, which have not been certified by Sartorius Stedim Systems.</p>	
Authorized person for documentation	Sartorius Stedim Biotek GmbH attn: Marc Höpner August-Spindler-Strasse 11 37079 Göttingen, Deutschland Phone +49 551 3083762, Fax +49 551 3082012	
Designation of the device	Biosafe freeze-thaw System	
Model version	S2 System, CryWedge System	
Cat.-No.	PTH-CS0000-0000, PTH-CS0000-0004, PTH-CS0000-0006, PTH-CS0000-0008, PTH-CS0000-0010, PTH-CS0000-0012, PTH-CS0000-0014, PTH-CS0000-0016	
Relevant directives of the EC	2004/42/EC Machinery 2004/108/EC Electromagnetic Compatibility 2006/95/EC Electrical equipment designed for use within certain voltage limits	
Applied harmonized standards	EN ISO 12100-1:2000, EN ISO 12100-2:2000 EN 61326-1:2006 EN 61010-1:2001	
Applied national standards and Technical specifications	not applied	
Date and Signature	24.01.2011	
Function of Signatory	Marc Höpner Director of R&D for Automation, Sensors and Instruments	 Dr. Susanne Dreihöfen Director of Quality Engineering Systems and Instruments

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Format:

```
CryoPilot 5.0\n
METADATA\n
Title\t<title>\n
Process Start\t<timestamp format=yyyy-MM-dd HH:mm:ssZ>\n
Computer\t<computer name>\n
User\t<user name>\n
Description\t<description>\n
START\n
Date/Time(UTC)\tElapsed Time(h)\tChiller Status\tHTF Supply Temperature Setpoint(C)\tHTF Supply Temperature(C)\tMixer Status\tProduct Temperature 1(C)\tProduct Temperature 2(C)\t...Product Temperature n(C)\n
<timestamp format=yyyy-MM-dd HH:mm:ssZ>\t<elapsedhours format=0.000000>\t<chiller = Off|On|Unknown>\t<value format=0.0>\t<value format=0.0>\t<mixer = Off|On|Unknown>\t<value format= 0.0>\t<value format= 0.0>\t...<value format=0.0>\n
.\n
.\n
.\n
.\n
.
```

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START

Date/Time(UTC) Elapsed Time(h) Chiller Status HTF Supply Temperature

Setpoint(C) HTF Supply Temperature(C) Mixer Status Product 1

Temperature(C) Product 2 Temperature(C) Product 3 Temperature(C)

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Copyright by Sartorius Stedim Biotech GmbH, Goettingen, Germany. All rights reserved. No part of this publication may be reprinted or translated in any form or by any means without the prior written permission of Sartorius Stedim Biotech GmbH. The status of the information, specifications and illustrations in this manual is indicated by the date given below. Sartorius Stedim Biotech GmbH reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Status:
August 2011,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany

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Sartorius Stedim Biotech S.A.