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
# Table of Contents

<b>Sterisart® NF gamma – Overview</b> .....	4
Product description .....	5
Packaging .....	5
Intended application .....	6
Specification .....	6
<b>1. Preparation for use in the sterilization area of an isolator</b> .....	7
<b>2. Perform these steps before filtration of the samples</b> .....	8
<b>3. Dissolving   dilution step</b> .....	8
<b>4. Sample filtration   washing the samples and filling the container with liquid culture media</b> .....	10

# Sterisart® NF gamma – Overview

Instructions for dissolving powders with Sterisart® NF 16470 to prepare samples for subsequent sterility test according to international pharmacopoeias

Pos.	Designation
1	Double-needle metal spike (3.7 cm) with yellow color coding
2	Double-needle metal spike (2.3 cm) without color coding
3	Double tubing, split in the middle to approx. 20 cm 1 tube with black plastic rings for identification
4	Additional plastic bag

### Scope of the Instructions for Use

These instructions are applicable to the following products:

#### 16470 GBD

Sterisart® NF gamma, tubing set for dissolving and diluting slightly soluble powders in closed containers (box of ten, double packed, sterilized by gamma irradiation)

Product description  
(see page 4)

The tubing system consists of two double-needle metal spikes (one 3.7 cm in length with yellow color coding and the other 2.3 cm in length) connected by double tubing. The double tubing is split half way, and one tube is identifiable by its black plastic rings.

Please note: The tube with the black rings does not go in the pump.

Packaging

#### **Sterisart® NF gamma:**

An individual package consists of a set of tubing in a gas-impermeable bag. This packaging allows the outer surfaces of the bag to be sterilized in the sterilization area of an isolator while preventing any sterilant (in this case H<sub>2</sub>O<sub>2</sub>) from penetrating through the packaging into the Sterisart® system. The label on the bag provides identifying information about the system (brief description, order number, expiration date, sterilization date, lot number).

A standard box contains ten sterile, single packed systems, which are additionally packed in a plastic bag. Each system has a sterilization indicator. The product has been reliably sterilized by gamma irradiation when the sterilization indicator changes from orange to red.

## Intended application

Sterisart® NF gamma systems have been specifically designed for sterility testing using the membrane filter method described in pharmacopoeias, in combination with a dedicated peristaltic pump. Type 16470 in combination with Sterisart® 16466 (4 cm double-needle metal spike for unvented containers) or type 16467 (6-cm metal needle) allows sterility testing of slightly soluble powders in closed, unvented containers.

Sterisart® NF gamma has been optimized for sterility testing in isolators.

## Important notes on safe and correct usage

**Sterisart® NF systems are only to be used with double-tubing peristaltic pumps designed for sterility testing with such systems, such as the Sterisart® Universal pumps 16419 and 16420.**

**All necessary precautions must be taken to avoid secondary contamination. Handle systems parts only when wearing sterile gloves.**

**The adapter needles are very sharp. Take great caution to avoid piercing gloves or injuring hands when handling them.**

**Check the color indicator and the expiration date on the box before removing a Sterisart® NF system from its packaging.**

## Specifications


Max. operating pressure:	3 bar at 20°C
Min. operating temperature:	50°C
Chemical compatibility:	Details available on request

Please follow the instructions for use starting on the next page.


# 1. Preparation for use in the sterilization area of an isolator

## **Sterisart® NF gamma**


Open the standard box and the plastic bag to remove the quantity of individual Sterisart® NF packages required for the number of sterility tests to be performed.

Place the individual packages of Sterisart® NF systems, still enclosed in gas-impermeable plastic bags, in the sterilization area.


Sterilize the bag according to the isolator validation specifications, for example, with H<sub>2</sub>O<sub>2</sub>.

Cut open the gas-impermeable bag to remove the set of tubing for the dissolving | dilution step.

## 2. Perform these steps before filtration of the samples


Cut open the bag as illustrated in 1.4 and remove the set of tubing.


Thread the tube without the black rings (not color coded) in the pump head, following the instructions given for the pump used.

Make sure that the yellow color-coded double-needle metal spike is located on the inlet side of the pump.


## 3. Dissolving | dilution step

Insert the double-needle metal spike as far as necessary into the closed container with solvent.


Then place the container upside down on the stand.


Insert the double-needle metal spike without color coding (the shorter of the two double-needle spikes) as far as necessary into the upright bottle filled with slightly soluble powder.

Switch on the pump to transfer the solvent into the container filled with powder (not illustrated). Thread the unmarked part of the tube only in the pump head, NOT the tube with the black rings.


The liquid flows into the bottle and the powder starts to dissolve slowly.

Once the bottle containing the powder is filled with solvent, the solution with the dissolved powder starts to flow back into the storage bottle.

In this manner, the solvent and the powder to be dissolved keep recirculating until the powder has dissolved completely.


After the powder has dissolved completely, turn both containers 180° and press the liquid out of the smaller one of the two bottles by pumping in air.

The dissolving step is completed once all the liquid has flowed back into the storage bottle (usually the larger of the two).

Dissolve additional samples from other sample containers according to the same procedure.

## 4. Sample filtration | washing the samples and filling the container with liquid culture media

Remove the double-needle metal spike of tubing set 16470 from the seal on the larger storage bottle.

Then filter the liquid in the storage bottle (with the dissolved powder) through a Sterisart® NF 16466 and | or 16467 system (the version has a spike with a vent filter). The procedure is detailed in the Instructions for Use for 16466 and | or 16467 alpha or gamma (Publication no. SLD6005- a).

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Copyright by  
Sartorius Stedim Biotech GmbH,  
Goettingen, Germany.  
All rights reserved. No part of  
this publication may be reprinted  
or translated in any form or  
by any means without the prior  
written permission of Sartorius  
Stedim Biotech GmbH.  
The status of the information,  
specifications and illustrations  
in this manual is indicated by  
the date given below.  
Sartorius Stedim Biotech GmbH  
reserves the right to make changes  
to the technology, features,  
specifications and design of the  
equipment without notice.

Status:  
August 2009,  
Sartorius Stedim Biotech GmbH,  
Goettingen, Germany