



Sartobind<sup>®</sup>  
Membrane Adsorbers for  
Rapid Purification of Proteins





# ... for Robust Separations

## Unique microporous structure

Sartobind Membrane Adsorbers display a macroporous structure with a pore size of  $> 3 \mu\text{m}$ .

That's orders of magnitudes larger than conventional chromatographic gel matrices. Molecules are transported by convective flow to ligands.

## Characteristics of Membrane Adsorbers (MA)

- Ready-to-use units
- Simple handling
- Pore sizes  $> 3$  and  $0.45 \mu\text{m}$
- Negligible diffusion limitation
- The disposable capsules with 4 mm bed height fit for polishing, the capsules with 8 mm bed height are for capturing and polishing (single-use | batch reuse). Both capsule groups are fully scalable
- Robust high performance separations
- No bed cracking, channeling, air entrapment
- Flow rate of ion exchange membranes 5–30 membrane volumes per minute
- Chemistries: strong and weak ion exchange, coupling, affinity and metal chelate ligands

## Low unspecific adsorption

The basis for all Sartobind membranes is a stabilized reinforced cellulose. It is made from regenerated cellulose and during the production to Sartobind it runs through a number of stabilization and grafting steps until a chromatographic matrix is formed on the cellulose backbone.

## Speed up

Sartobind ready-to-use units are run at 5 to 30 membrane volumes per minute. This is at least one order of magnitude faster than chromatographic columns. A typical speed up factor is about 25 measured in direct comparison to conventional column technology (reference: Walter, J. K. in: *Bioseparation and Bioprocessing, Strategies and Considerations for Advanced Economy in Downstream Processing of Biopharmaceutical Proteins*. G. Subramanian (ed.) Wiley VCH, Vol. II p. 447–460 (1998). Flow rate does not affect break through performance (see fig. 2).




Fig. 2: Break through performance of a Sartobind 35 mL cation exchange module at monoclonal antibody purification with various flow rates. For other published data see also: H.L. Knudsen et al., Genentech Inc., *J. Chromatogr. A* 907 (2001) 145–154.

## Sartobind membrane types

- Sartobind S, Q, STIC PA and D ion exchangers
- Sartobind IDA (iminodiacetic acid) metal chelate
- Sartobind Aldehyde
- Sartobind Protein A (recombinant)

## Sartobind applications

### Purification and Concentration

- Viruses, large proteins, blood factors, protein conjugates, VLP, mAbs

### Contaminant removal

- Host cell proteins, DNA, viruses, leached ligands, endotoxins, aggregates

## Chemical compatibility

The housing of Sartobind 4 and 8 mm capsules is polypropylene. The housing of Sartobind MA 15, 75 and 100 is polysulfone which is stable to many standard solvents applied in chromatography. Sartobind ion exchange membranes are compatible with alcohols such as ethanol, isopropanol, glycerol, and denaturing solvents such as 8 M urea and 8 M guanidine HCl and can be cleaned with 1 N sodium hydroxide.

## Constant capacity

The robustness of Sartobind Membrane Adsorbers in ion exchange and affinity chromatography has been tested in consecutive runs of 1000 chromatographic cycles (Fig. 3). Sartobind can be reused many times, and affinity membranes maintain their binding capacity over many cycles reflecting the chemical stability of the ligands.




Fig. 3: Purification cycles of bovine serum albumin on Sartobind Q 75 strong anion exchanger repeated 1000 times. Flow rate: 10 mL/min, cycle time: 10 min, equilibration buffer: 20 mM phosphate buffer pH 7.0, sample: 5 mL bovine serum diluted 1:20 with equilibration buffer, elution buffer: 20 mM phosphate buffer pH 7.0 + 1 N NaCl, regeneration after each 100 cycles with 1 N NaOH

## Validated Technology




Fig. 4: There are no restrictions regarding ligands and applications.

### Quality assurance

Membranes are tested for protein dynamic binding capacity, charge density, thickness and flow rate. Capsules are tested for dynamic binding capacity and flow rate.

### Validated technology

Sartobind 4 mm and 8 mm capsules have met the requirements of current USP Plastic Class V and VI, respectively. Validation guides are available with detailed technical data and test results according to USP. Additionally, extractables guide is available for modules with data of GC-MS, HPLC and FT-IR analysis from aqueous, ethanolic and sodium hydroxide extraction.

### FDA approved

Sartobind is the first Membrane Adsorber applied in approved pharmaceutical production (approval of CAMPATH-1H<sup>®</sup> by Food and Drug Administration, March 2001).

## Formats from Lab to Process



Fig. 5: Sartobind MAs may be used by hand or with a chromatographic system via Luer adapters.



Fig. 6: Sartobind MA units for laboratory applications



Fig. 7: Sartobind Protein A 2 mL runs at 10 mL/min with a dynamic binding of 10–15 mg IgG per unit

### Sartobind MA units

You may use Sartobind Membrane Adsorbers for any ion exchange or affinity chromatography which require high speed and simple operation. The small scale units for laboratory use is the MA 15 with 0.8 mm bed height (3 membrane layers), MA 75 with 4 mm bed height (15 layers) and MA 100 with 1.4 mm (5 layers). These may be reused hundreds of times.

### Work horse Sartobind MA 75 for lab purification work

The Sartobind MA 75 format with 4 mm bed height is available with all membrane types. It can be scaled straightly to the production format with the same bed height. They are designed for high throughput laboratory purification of biomolecules in the range of 5–50 mg.

### Protein A

The Sartobind Protein A is the ideal tool for ultra rapid quantification of immunoglobulins (IgG). You may achieve much more sample throughput per mL than with any other Protein A matrix. You can use the unit by a syringe and perform the quantification within 10 minutes. The enclosed Luer adapters allow direct connection to a liquid chromatography system.



Fig. 9: Sartobind pico

### Sartobind pico 0.08 mL

This format keeps the bed height of larger 4 mm capsules and allows for small scale trials with minute amount of samples e.g. for virus removal studies.



Fig. 10: Vivapure centrifugal units

### Centrifugal format

If you don't need flow control but you have to screen a large number of proteins in parallel, use the Vivapure<sup>®</sup> centrifugal units.



Fig. 11: 8 strip | 96 well plates

### 8 strips in 96 well plates

This format is suitable for process developers in biopharma to screen a large number of buffer conditions in a short time.



Fig. 8: Sartobind capsules with 8 mm bed height: for single-use or batch-use bind and elute or flow through applications

### **Production: Sartobind Jumbo reusable and disposable**

The Jumbo is a 5 liter membrane chromatography device for large scale capturing and impurity removal at high flow rates. The Sartobind nano 3 mL is the appropriate scale down model with same bed height of 8 mm.



Fig. 12: Sartobind 150 mL capsule cross section displays small upstream and downstream channels

### **Void volume optimized construction**

Sartobind 8 mm capsules are designed to achieve smallest peak | pool volumes. That is the reason for sharp peak resolution of this capsule family and the applicability for bind and elute purification processes.



## Sartobind 4 mm Capsules for Contaminant Removal




Fig. 13: Sartobind 4 mm capsule scale up from mL to 10000 liter scale: Sartobind pico 0.08 mL to Sartobind mega 1.62 L

### High speed polishing

Sartobind 4 mm capsules are designed to remove charged contaminants from therapeutic proteins at accelerated flow rates of 10 to 30 membrane volumes per minute by ion exchange membrane chromatography. The high throughput is a direct result of negligible mass transfer effects and is made possible by the  $>3 \mu\text{m}$  macroporous membrane with 4 mm bed height.

### Capsule design

The Sartobind 4 mm looks like a standard filter capsule except that the adsorbing membrane is rolled up on a core to form a cylinder. The flow is from top into the outside channel and then perpendicular through the membrane layers to the center of the cylinder and leaves the capsule at the outlet.

### Downscale with Sartobind pico and nano

The Sartobind pico is the smallest member of the 4 mm capsule family with 0.08 mL membrane volume scalable up to the largest capsules. It saves costs in sample consumption and virus spiking studies. Sartobind nano 1 mL is the smallest cylindric capsule with 4 mm bed, Sartobind pico displays a flat sheet design to obtain the smallest scalable device.

### Membrane ligands available with 4 mm capsules

- Quaternary ammonium - anion exchanger
- Sartobind STIC primary amine - salt tolerant anion
- Sulfonic acid - cation exchanger

### Efficient

- Higher throughput for trace impurity removal at 10 kg protein per liter membrane in flow through mode

### Economical

- No hardware investment & maintenance
- No column packing, testing, regeneration
- No re-use validation
- Less labor and smaller footprint
- Buffer consumption may be decreased 95% compared to columns
- Sartobind STIC PA reduces tank investments and buffer usage as cation exchange pools are cleared from contaminants without dilution

### Easy to use

- Disposable
- Simple and fast set up
- Handling like a filter capsule
- Less complex

## A Standard for Virus Purification

### Unique for large proteins, DNA and viruses

Sartobind membranes are unique for purification of large proteins and virus particles. Conventional gel beads possess 95% of the binding capacity within the bead and large particles cannot enter their small pores. In contrast Sartobind membranes display huge capacities for large molecules. E.g. for DNA 10 times more dynamic capacity than conventional beads (reference: Sartobind Application Note DNA removal) and for virus purifications gels cannot compete at all. Sartobind membranes have been successfully applied for the purification of alphaherpesvirus, adenovirus and adenoassociated virus. The Vivapure AdenoPACK adenovirus purification and concentration kits offer researchers a convenient way to clean up within 1-2 hours. The LentiSELECT kit reduces the purification time of lenti viruses up to  $8 \times 10^6$  particles below 1 hour. Both systems are scalable to production format.




Fig. 14: Viruses can bind to the internal membrane structure



Fig. 15: Adeno PACK for convenient adenovirus purification and concentration

# Technical Data

To connect to a chromatography system or Luer syringe



	pico	MA 15	MA 75   2mL	MA 100	nano
Application	Down-scale for biopharma development	Screening, application optimization, reusable	Screening, application optimization, reusable	Screening, application optimization, reusable	Down-scale for biopharma production
Membrane material	Stabilized reinforced cellulose				
Adsorption area [cm <sup>2</sup> ]	2.9	15	75   100	100	36   110
Number of layers	15	3	15   20	5	15   30
Bed height [mm]	4	0.8	4.0	1.4	4   8
Bed volume* [mL]	0.08	0.41	2.1   2	2.8	1   3
Membrane diameter [mm]	5	25	25	50	–
Housing material	Polypropylene	Polysulfone	Polysulfone	Polysulfone	Polypropylene
Inlet connector	Female Luer	Female Luer	Female Luer	Female Luer	Female Luer
Outlet connector	Female Luer	Male Luer	Male Luer	Male Luer	Female Luer
Typical dynamic binding capacity 10% binding capacity* [mg/unit]	2.3 / Q pico 4.0 / PA pico 2.0 / S pico 0.6 / Phenyl pico	10.5 / S 15 12 / Q 15	52.5 / S 75 60 / Q 75 60 / D 75 7.5 / IDA 75 15 / Protein A	70 / S 100 80 / Q 100	29 / Q nano 1 mL 25 / S nano 1 mL 60 / PA nano 1 mL 88 / Q nano 3 mL 77 / S nano 3 mL 44 / Phenyl nano 3 mL
Recommended flow rate [mL/min]	0.8–2.4	10–30	10–30	20–50	10–30   5–15
Maximum pressure [MPa]	0.6	0.6	0.6	0.6	0.4

\* Reference proteins ion exchange: S: lysozyme, Q, D, PA: BSA, IDA: 6 × His-tagged protein, epoxy. IgG, Protein A: IgG, Phenyl: IgG

## Overview Membrane Types and Formats


Membrane Type	Description	Pore size (µm)	Capacity* (mg/mL)	Sartobind Lab Scale					Process Scale	
				A4 sheet	96 well 19 µl	MA 15, 100	MA 75	pico 0.08 mL	Capsule FT 4 mm	Capsule FT/B&E** 8 mm
S	Sulfonic acid	>3	25	•	•	•	•	•	•	•
Q	Quaternary ammonium	>3	29	•	•	•	•	•	•	•
PA	Primary amine	>3	60		•			•	•	
D	Diethylamine	>3	22	•			•			
IDA	Iminodiacetic acid	>3	3.6	•			•			
HIC	Phenyl	>3	44		•			•		•
Aldehyde	Aldehyde	0.45	1	•						
Protein A	Recombinant Protein A	0.45	5–7.5				•			

\* Typical dynamic binding capacity 10%

\*\* Flow through | Bind and elute

• available

**Spin colmuns for operation with centrifuge**



**Mini**

**Maxi**

Application	Screening, application optimization, disposable	Screening, application optimization, disposable
Membrane ligands	Q, D, S	Q, S
Membrane diameter [mm]	6.4	22.5
Housing material	Polypropylene	Polypropylene
Minimum static binding capacity* [mg/unit]	H: 4	H: 48–64
Maximum volume [mL] (device inner volume)	H: 0.4	H: 19
Spin speed [g]	H: 2000	H: 500
Storage before use at	Room temperature	Room temperature

M: medium, H: high

\* Reference proteins ion exchange: S: lysozyme, Q, D: BSA

**Screening devices for operation with centrifuge | vacuum**



**96-well plates equipped with 8 strips**

Application	Screening for process development work
Membrane ligands	S, Q, STIC PA and HIC Phenyl
Housing material	Polypropylen
Number of layers	3
Bed height	0.8 mm
Membrane volume	0.019 mL
Membrane diameter	5.5 mm
Maximum well volume [mL]	0.5 mL
Collection plate well volume [mL]	2 mL

# Ordering Information

## Luer Lock units | sheets | kit

Order No.	Sartobind MA units   Kit	Quantity
93IEXQ42GB-12--A	Sartobind Q 15	4
93IEXS42GB-12--A	Sartobind S 15	4
93IEXQ42DB-12--V	Sartobind Q 75	2
93IEXS42DB-12--V	Sartobind S 75	2
93IEXD42DB-12--V	Sartobind D 75	2
93IEXQ42BC-12	Sartobind Q 100	1
93IEXS42BC-12	Sartobind S 100	1
93IDA426DB-12--V	Sartobind IDA 75	2
93PRAP06HB-12--A	Sartobind Protein A 2 mL	4
90-KIT-01	Sartobind Demo Kit	

Order No.	Sartobind A4 Sheet	Quantity
94IEXS42-001	Sartobind S A4 Sheet	1
94IEXQ42-001	Sartobind Q A4 Sheet	1
94IEXD42-001	Sartobind D A4 Sheet	1
94IDA-42-001	Sartobind IDA A4 Sheet	1
94ALD-06-001	Sartobind Aldehyde A4 Sheet	1

## Sartobind capsules 4 mm bed height

Order No.	Sartobind Q	Quantity
92IEXQ42DD-11--D	Sartobind Q pico 0.08 mL	10
92IEXQ42DN-11	Sartobind Q SingleSep nano 1 mL	1
92IEXQ42DN-11--A	Sartobind Q SingleSep nano 1 mL	4
92IEXQ42D4-00--A	Sartobind Q SingleSep mini capsule, hose barb	4
92IEXQ42D4-SS--A	Sartobind Q SingleSep mini capsule, sanitary	4
92IEXQ42D9-00--A	Sartobind Q SingleSep 5" capsule, hose barb	4
92IEXQ42D9-SS--A	Sartobind Q SingleSep 5" capsule, sanitary	4
92IEXQ42D1-SS	Sartobind Q SingleSep 10"capsule	1
92IEXQ42D2-SS	Sartobind Q SingleSep 20"capsule	1
92IEXQ42D3-SS	Sartobind Q SingleSep 30"capsule	1
92IEXQ42DC3SS	Sartobind Q SingleSep mega capsule	1

Order No.	Sartobind S	Quantity
92IEXS42DD-11--D	Sartobind S pico 0.08 mL	10
92IEXS42DN-11	Sartobind S SingleSep nano 1 mL	1
92IEXS42DN-11--A	Sartobind S SingleSep nano 1 mL	4
92IEXS42D4-00--A	Sartobind S SingleSep mini capsule, hose barb	4
92IEXS42D4-SS--A	Sartobind S SingleSep mini capsule, sanitary	4
92IEXS42D9-00--A	Sartobind S SingleSep 5" capsule, hose barb	4
92IEXS42D9-SS--A	Sartobind S SingleSep 5" capsule, sanitary	4
92IEXS42D1-SS	Sartobind S SingleSep 10"capsule	1
92IEXS42D3-SS	Sartobind S SingleSep 30"capsule	1

Order No.	Sartobind STIC	Quantity
92STPA42DD-11--D	Sartobind STIC PA pico 0.08 mL	10
92STPA42DN-11--A	Sartobind STIC PA nano 1 mL	4
92STPA42D9-FF--A	Sartobind STIC PA 5" capsule	4
92STPA42D1-SS	Sartobind STIC PA 10" capsule	1
92STPA42D2-SS	Sartobind STIC PA 20"capsule	1
92STPA42D3-SS	Sartobind STIC PA 30" capsule	1
92STPA42DC3SS	Sartobind STIC PA mega capsule	1

## Sartobind capsules 8 mm bed height

Order No.	Sartobind Q	Quantity
96IEXQ42EUC11--A	Sartobind Q nano 3 mL	4
96IEXQ42E9BFF	Sartobind Q 150 mL	1
96IEXQ42E3FSS	Sartobind Q 1.2 L	1
96IEXQ42E3ESS	Sartobind Q Jumbo 5 L	1

Order No.	Sartobind S	Quantity
96IEXS42EUC11--A	Sartobind S nano 3 mL	4
96IEXS42E9BFF	Sartobind S 150 mL	1
96IEXS42E3FSS	Sartobind S 1.2 L	1
96IEXS42E3ESS	Sartobind S Jumbo 5 L	1

Order No.	Sartobind Phenyl	Quantity
92HICP42DD-11--D	Sartobind Phenyl pico 0.08 mL	10
96HICP42EUC11--A	Sartobind Phenyl nano 3 mL	4
96HICP42E9BFF	Sartobind Phenyl 150 mL	1
96HICP42E3FSS	Sartobind Phenyl 1.2 L	1
96HICP42E3ESS	Sartobind Phenyl Jumbo 5 L	1

#### Virus purification kits

Order No.	AdenoPACK   LentiSELECT	Quantity
VS-AVPQ020	Vivapure AdenoPack 20	
VS-AVPQ022	Vivapure AdenoPack 20 RT	
VS-AVPQ101	Vivapure Adenopack 100	
VS-AVPQ102	Vivapure Adeno PACK 100 RT*	
VS-AVPQ501	Vivapure Adeno PACK 500	
VS-AVPQ502	Vivapure AdenoPack 500 RT	
VS-LVPQ040	Vivapure LentiSELECT 40	

\* Vivapure Adeno PACK RT does not contain Benzonase®

#### Accessories

Order No.	Adapters   Holders	Quantity
1ZA---0004	Adapter UNF 10-32 to Luer male, PEEK	1
1ZA---0005	Adapter UNF 10-32 to Luer female, PEEK	1
17002---140	Pair of Luer adapters, black Tefzel M6 thread	2
16517-----E	Syringe filter holders, 25 mm diameter, polycarbonate for filter or adsorber membrane	12
16214	Syringe filter holder 25 mm diameter, stainless steel for filter or adsorber membrane	1
9ZAIAM0001	Legs for Sartobind mega capsule (stainless steel)	3
9ZGL--0102	Trolley for Sartobind Jumbo 5 L capsule	1

#### Spin columns

Order No.	Vivapure Mini Spin Columns	Spin columns   Centrifuge tubes
VS-IX01SH24	Vivapure S Mini H	24   48
VS-IX01QH24	Vivapure Q Mini H	24   48
VS-IX01DH24	Vivapure D Mini H	24   48

Order No.	Vivapure Maxi Spin Columns	Spin columns   Centrifuge tubes
VS-IX20SH08	Vivapure S Maxi H	8   16
VS-IX20QH08	Vivapure Q Maxi H	8   16

#### 96-well plate

Order No.	96-well plates   Accessories   Software	Quantity
99IEXQ42GC-----V	Sartobind® Q 96-well plate, 2 units	2 (24 × 8-strips)
99IEXQ42GC-----D	Sartobind® Q 96-well plate, 10 units	10 (120 × 8-strips)
99STPA42GC-----V	Sartobind® STIC PA 96-well plate, 2 units	2 (24 × 8-strips)
99STPA42GC-----D	Sartobind® STIC PA 96-well plate, 10 units	10 (120 × 8-strips)
99IEXS42GC-----V	Sartobind® S 96-well plate, 2 units	2 (24 × 8-strips)
99IEXS42GC-----D	Sartobind® S 96-well plate, 10 units	10 (120 × 8-strips)
99HICP42GC-----V	Sartobind® Phenyl 96-well plate, 2 units	2 (24 × 8-strips)
99HICP42GC-----D	Sartobind® Phenyl 96-well plate, 10 units	10 (120 × 8-strips)
VW96VAC01	Vac96 vacuum manifold	1
VW96VAA02	Vac96 liquid trap and reservoir	1
VW96VAA04	96 deep well collection plate 2 mL (square wells)	25
VW96VAA05	Replacement seal for Vac96 vacuum manifold	25
VW08VAC01	Vac8 vacuum manifold	1
VW08VAA02	Vac8 liquid trap and reservoir	1
VW08VAA03	8-well collection strips 1.2 mL (round wells)	125
VW08VAA04	Replacement seal for Vac8 vacuum manifold	1
BBI-8582670	Swing-out rotor for deep well plates (max. 64 mm high) incl. 2 buckets	1
16694-2-50-06	Microsart® mini.vac vacuum pump 230 V   50 Hz	1
16694-1-60-06	Microsart® mini.vac vacuum pump 115 V   60 Hz	1
BB-8870600	MODDE Single User License (with 1 year maintenance)	1
BB-8870602	MODDE Single University License (with 1 year maintenance, 85% discount)	1

# Sales and Service Contacts

For further contacts, visit [www.sartorius-stedim.com](http://www.sartorius-stedim.com)

## Europe

### Germany

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen

Phone +49.551.308.0  
Fax +49.551.308.3289

Sartorius Stedim Systems GmbH  
Robert-Bosch-Strasse 5 – 7  
34302 Guxhagen

Phone +49.5665.407.0  
Fax +49.5665.407.2200

### France

Sartorius Stedim FMT S.A.S.  
ZI Les Paluds  
Avenue de Jouques – CS 91051  
13781 Aubagne Cedex

Phone +33.442.845600  
Fax +33.442.845619

Sartorius Stedim France SAS  
ZI Les Paluds  
Avenue de Jouques – CS 71058  
13781 Aubagne Cedex

Phone +33.442.845600  
Fax +33.442.846545

### Austria

Sartorius Stedim Austria GmbH  
Franzosengraben 12  
1030 Vienna

Phone +43.1.7965763.18  
Fax +43.1.796576344

### Belgium

Sartorius Stedim Belgium N.V.  
Leuvensesteenweg, 248/B  
1800 Vilvoorde

Phone +32.2.756.06.80  
Fax +32.2.756.06.81

### Hungary

Sartorius Stedim Hungária Kft.  
Kagyló u. 5  
2092 Budakeszi

Phone +36.23.457.227  
Fax +36.23.457.147

### Italy

Sartorius Stedim Italy S.p.A.  
Via dell'Antella, 76/A  
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41  
Fax +39.055.63.40.526

### Poland

Sartorius Stedim Poland Sp. z o.o.  
ul. Wrzesinska 70  
62-025 Kostrzyn

Phone +48.61.647.38.40  
Fax +48.61.879.25.04

### Russian Federation

LLC "Sartorius ICR"  
Uralskaya str. 4, Lit. B  
199155, Saint-Petersburg

Phone +7.812.327.5.327  
Fax +7.812.327.5.323

### Spain

Sartorius Stedim Spain SA  
C/Isabel Colbrand 10,  
Oficina 70  
Polígono Industrial de Fuencarral  
28050 Madrid

Phone +34.90.2110935  
Fax +34.91.3589623

### Switzerland

Sartorius Stedim Switzerland AG  
Ringstrasse 24 a  
8317 Tagelswangen

Phone +41.52.354.36.36  
Fax +41.52.354.36.46

### U.K.

Sartorius Stedim UK Ltd.  
Longmead Business Centre  
Blenheim Road, Epsom  
Surrey KT19 9 QQ

Phone +44.1372.737159  
Fax +44.1372.726171

## America

### USA

Sartorius Stedim North America Inc.  
5 Orville Drive, Suite 200  
Bohemia, NY 11716

Toll-Free +1.800.368.7178  
Fax +1.631.254.4253

### Argentina

Sartorius Argentina S.A.  
Int. A. Ávalos 4251  
B1605ECS Munro  
Buenos Aires

Phone +54.11.4721.0505  
Fax +54.11.4762.2333

### Brazil

Sartorius do Brasil Ltda  
Av. Dom Pedro I, 241  
Bairro Vila Pires  
Santo André  
São Paulo  
Cep 09110-001

Phone +55.11.4451.6226  
Fax +55.11.4451.4369

### Mexico

Sartorius de México S.A. de C.V.  
Circuito Circunvalación Poniente  
No. 149  
Ciudad Satélite  
53100, Estado de México  
México

Phone +52.5555.62.1102  
Fax +52.5555.62.2942

## Asia | Pacific

### Australia

Sartorius Stedim Australia Pty. Ltd.  
Unit 5, 7-11 Rodeo Drive  
Dandenong South Vic 3175

Phone +61.3.8762.1800  
Fax +61.3.8762.1828

### China

Sartorius Stedim Biotech (Beijing) Co. Ltd.  
No. 33 Yu'an Road  
Airport Industrial Park Zone B  
Shunyi District, Beijing 101300

Phone +86.10.80426516  
Fax +86.10.80426580

Sartorius Stedim Biotech (Beijing) Co. Ltd.  
Shanghai Branch Office  
3rd Floor, North Wing, Tower 1  
No. 4560 Jin Ke Road  
Pudong District, Shanghai 201210

Phone +86.21.68782300  
Fax +86.21.68782332 | 68782882

Sartorius Stedim Biotech (Beijing) Co. Ltd.  
Guangzhou Representative Office  
Unit K, Building 23  
Huihua Commerce & Trade Building  
No. 80 Xianlie Middle Road  
Guangzhou 510070

Phone +86.20.37618687 | 37618651  
Fax +86.20.37619051

### India

Sartorius Stedim India Pvt. Ltd.  
#69/2-69/3, NH 48, Jakkasandra  
Nelamangala Tq  
562 123 Bangalore, India

Phone +91.80.4350.5250  
Fax +91.80.4350.5253

### Japan

Sartorius Stedim Japan K.K.  
4th Fl., Daiwa Shinagawa North Bldg.  
8-11, Kita-Shinagawa 1-chome  
Shinagawa-ku, Tokyo, 140-0001 Japan

Phone +81.3.4331.4300  
Fax +81.3.4331.4301

### Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.  
Lot L3-E-3B, Enterprise 4  
Technology Park Malaysia  
Bukit Jalil  
57000 Kuala Lumpur, Malaysia

Phone +60.3.8996.0622  
Fax +60.3.8996.0755

### Singapore

Sartorius Stedim Singapore Pte. Ltd.  
1 Science Park Road,  
The Capricorn, #05-08A,  
Singapore Science Park II  
Singapore 117528

Phone +65.6872.3966  
Fax +65.6778.2494

### South Korea

Sartorius Korea Biotech Co., Ltd.  
8th Floor, Solid Space B/D,  
PanGyoYeok-Ro 220, Bundang-Gu  
SeongNam-Si, GyeongGi-Do, 463-400

Phone +82.31.622.5700  
Fax +82.31.622.5799



◀ [www.sartorius-stedim.com](http://www.sartorius-stedim.com)